

Pathogen ID & Surveillance: Sequencing Roadmap Beyond COVID-19 at illumina

Let's Discuss!!!

Chris Detter — cdetter@illumina.com (505-699-1996)

AMR East & LATAM; Microbial and Infectious Disease Sales Specialist

Genome Trakr Meeting – Technology Talk

College Park, MD

10.20.2022

For Research Use Only. Not for use in diagnostic procedures.

illumina[®]

© 2021 Illumina, Inc. All rights reserved.



illumina Infectious Disease & Microbiology Specialist Team

- **Chris Detter** (505-699-1996; cdetter@illumina.com)

ID&M Specialist – US East and LATAM

- DoE Joint Genome Institute & Livermore National Lab
 - Part of Human Genome Program... R&D
- Los Alamos National Lab – Genome Center Director
 - Focus on Infectious Disease Genomics....and SFAF
- Defense Threat Reduction Agency (DTRA)
 - Built Genome centers abroad
- MRIGlobal
 - Chief Scientist R&D for Biosurveillance and Infectious Disease
 - Applied development of Biosurveillance using Genomics



- **Jason Smith** (858-531-8182; jsmith6@illumina.com)

ID&M Specialist – US West and Canada

- Illumina: Sequencing Specialist – US-West
- Past Companies: PacBio, Ion Torrent, Applied Bio, LLNL
- Location: Reno, NV
- Expertise: Microbiology (UCSB) , Computer Science (CSU Chico), Next-Generation-Sequencing (Field Sales/Support)



Emerging and Reemerging Infectious Diseases

List of priority infectious diseases that threaten global health security

Viral infections

Pandemic influenza

COVID-19¹

Middle East respiratory syndrome¹

Severe acute respiratory syndrome¹

Marburg¹

Ebola¹

Lassa fever¹

Viral hemorrhagic fevers

Rift Valley fever¹

West Nile fever

Dengue

Polio

Zika¹

Nipah and henipavirus diseases¹

Chikungunya¹

Measles

Yellow fever

Viral hepatitis

Bacterial infections

Tuberculosis¹

Invasive meningococcal disease

Invasive pneumococcal disease

Cholera

Typhoid

Diphtheria

Pertussis (whooping cough)

Drug-resistance

Drug-resistant bacterial, viral, and protozoal infections

Drug-resistant malaria


Antiretroviral-resistant human immunodeficiency virus/AIDS




1 - WHO Blueprint priority disease.

Zumla A, Hui DSC. Emerging and Reemerging Infectious Diseases: Global Overview. *Infect Dis Clin North Am*. 2019;33(4):xiii-xix. doi:10.1016/j.idc.2019.09.001

2022 Infectious Disease & Microbial Focus at Illumina



Infectious Disease Testing



Public Health



Microbiology Research


Illumina Developed

Partnerships

Commercial Enablement

Next-Generation Sequencing Unlocking the Power of the Microbial Genome

> APPLICATIONS



Illumina technology has been used in the fight against COVID-19 at all stages of the pandemic




JANUARY 2020

IDENTIFICATION OF NEW
VIRUS CAUSING COVID-
19 WITH ILLUMINA
PLATFORM



JUNE 2020

ILLUMINA DEVELOPS FIRST
SEQUENCING BASED
COVID-19 DIAGNOSTIC TEST
COVIDSeq™



SINCE LATE 2020

PUBLIC HEALTH AND
RESEARCH INSTITUTIONS USE
ILLUMINA TECHNOLOGY TO
CHARACTERIZE NOVEL
SARS-COV-2 VARIANTS
COVIDSeq, RVOP, RPIP, etc.

Discovery



Detection

Surveillance

Beyond Positive Detection of Sars-CoV-2 by RT-qPCR or NGS



Are there mutations in the viral sequence?

- Guidance for PCR-based testing revision
- Guidance for vaccine targets



What was the route of transmission?

- Guidance for public health policy decisions



Source: GISAID

How fast is the virus mutating?

Are there any co-infections?

[illegible]

<https://brcf.medicine.umich.edu/cores/advanced-genomics/faqs/sanger-sequencing-faqs/interpretation-of-sequencing-chromatograms/>. <https://arep.med.harvard.edu/Polonator/>

The Illumina RNA (DNA) Prep with Enrichment Workflow


Leveraging BLT technology to achieve rapid turnaround time without compromising quality

cDNA Synthesis



cDNA synthesis kit
included
(RNA)

Tagmentation



BLT rapid library prep
generation

Capture



3-Plex, rapid single-
hybridization

Quant & Norm



Completed libraries
ready to load on any
Illumina sequencer

Library Preparation*

9.5 (6.5)
hours

Comparable to Illumina DNA Prep w/
Enrichment + cDNA Synthesis

MiniSeq Rapid 100-cycles kit



<5h

1x100
18 samples

MiSeq V3 150-cycles kit



~12h 21h

1x75 2x75
24 samples

NextSeq 550 Mid-Output 150-cycles kit



~9h 15h

1x75 2x75
96 samples

For Research Use Only. Not for use in diagnostic procedures

*Approximate time for manual processing of 24
samples with 3-plex enrichment is ~9.5 from Total
RNA to library with ~2 hours of hands-on time.



RVOP

For Research Use Only. Not for use in diagnostic procedures

illumina®

1

Respiratory Viral Oligo Panel (RVOP)


Comprehensive Detection of 40+ Respiratory Viral Whole Genomes

- 7,800 probes (80bp) tiling the full genome of common respiratory viruses, recent flu strains, and SARS-CoV-2
- Includes probes that target human mRNAs to act as positive controls




Respiratory Viruses Detected by Respiratory Virus Oligo Panel


- Human coronavirus 229E
- Human coronavirus NL63
- Human coronavirus OC43
- Human coronavirus HKU1
- SARS-CoV-2
- Human adenovirus B1
- Human adenovirus C2
- Human adenovirus E4
- Human bocavirus 1 (Primate bocaparvovirus 1 isolate st2)
- Human bocavirus 2c PK isolate PK-5510
- Human bocavirus 3



- Human parainfluenza virus 1
- Human parainfluenza virus 2
- Human parainfluenza virus 3
- Human parainfluenza virus 4a
- Human metapneumovirus (CAN97-83)
- Respiratory syncytial virus (type A)
- Human Respiratory syncytial virus 9320 (type B)
- Influenza A virus (A/Puerto Rico/8/1934(H1N1))
- Influenza A virus (A/Korea/426/1968(H2N2))
- Influenza A virus (A/New York/392/2004(H3N2))
- Influenza A virus (A/goose/Guangdong/1/1996(H5N1))
- Human bocavirus 4 NI strain HBoV4-NI-385
- KI polyomavirus Stockholm 60
- WU Polyomavirus
- Human parechovirus type 1 PicoBank/HPeV1/a
- Human parechovirus 6
- Human rhinovirus A89
- Human rhinovirus C (strain 024)
- Human rhinovirus B14
- Human enterovirus C104 strain: AK11
- Human enterovirus C109 isolate NiCA08-4327



- Influenza A virus (A/Zhejiang/DJID-ZJU01/2013(H7N9))
- Influenza A virus (A/Hong Kong/1073/99(H9N2))
- Influenza A virus (A/Texas/50/2012(H3N2))
- Influenza A virus (A/Michigan/45/2015(H1N1))
- Influenza B virus (B/Lee/1940)
- Influenza B virus (B/Wisconsin/01/2010)
- Influenza B virus (B/Brisbane/60/2008)
- Influenza B virus (B/Colorado/06/2017)
- Influenza B virus (B/Washington/02/2019)
- Human control genes





RPIP

For Research Use Only. Not for use in diagnostic procedures

illumina®

2

Respiratory Pathogen ID/AMR Panel (RPIP)

Comprehensive Detection of Respiratory Pathogens and Antimicrobial Resistance Genes

Explore Complete List of Bacteria, Fungi, and Viruses on our webpage

Powered by



>280 Pathogens



180+
Bacteria



40+
Viruses*



50+
Fungi

>16 Antibiotic Classes




50+
Antimicrobial
Resistant Genes


1200+
Antimicrobial
Resistance Alleles

*full genome coverage available for SARS-CoV-2 and flu viruses

PDF Report



Computer Parsable Format (JSON)



NEW

VSP

For Research Use Only. Not for use in diagnostic procedures

illumina®

3

Introducing Viral Surveillance Panel (VSP)



- **66 Viral species** with 216 strains total (numbers in parenthesis)
- Oligo Hybrid Capture-80 bp oligos tiles across each genome for WGS Sequencing in **full-length**
- Flexible to mutations/strain **variation of 20%** due to hybrid capture technology
- Free to use DRAGEN Microbial Enrichment BaseSpace® app for analysis- including, consensus sequences, contig assembly, genome coverage and read counts
- Designed for **broad viral genomic surveillance for a variety of sample types*** including wastewater, environmental and post-clinical for public health



Arthropod transmitted / febrile tropical

Chikungunya(1)
Dengue Virus 1(1)
 Dengue Virus 2(1)
 Dengue Virus 3(1)
 Dengue Virus 4(1)
 West Nile Virus(2)
 Yellow fever virus(1)
 Zika(2)

Bloodborne

Hepatitis A(1)
 Hepatitis B (1)
 Hepatitis C(7)
 Hepatitis E(1)
 HIV 1(1)
 HIV 2(1)
 Torque Teno virus(1)

Encephalitis

Eastern equine encephalitis virus(1)
 Hendra
 henipavirus(1)
 Japanese encephalitis Virus(1)
 Nipah Virus(1)
 Tick-borne encephalitis virus(1)
 Venezuelan equine encephalitis virus(1)
 Western equine encephalitis virus(1)

Respiratory /Cardiopulmonary

Adenovirus(3)
 Coronavirus -229E(1)
 Coronavirus -HKU1(1)
 Coronavirus- OC43(1)
 Coronavirus-NL63(1)
 Influenza A(7)
 Influenza B(1)
 Hantavirus(33)
 MERS-CoV(1)
 Metapneumovirus(1)
 Parainfluenza(4)
 Parechovirus(1)
 Respiratory syncytial virus(1)
 Rhinovirus(3)
 SARS-COV(1)
SARS-COV-2(1)

Hemorrhagic fever

Chapare virus(1)
 Crimean-congo haemorrhagic fever virus(1)
Ebola(1)
 Guanarito Virus(1)
 Rift Valley Fever Virus(3)
 Junin Virus(1)
 Kyasanur Forest disease virus(1)
 Lassa fever virus(2)
 Lujo Hemorrhagic Fever Virus(1)
 Machupo Virus(1)
Marburg Virus(1)
 Omsk hemorrhagic fever virus(1)
 Sabia virus(1)

Enteric

Aichivirus (3)
 Astrovirus(1)
 Coxsackieviruses(28)
 Enterovirus(28)
 Norovirus(14)
Poliovirus(1)
 Rotavirus(4)
 Salivirus(1)
 Sapovirus(1)

Oncolytic

Oncolytic HPV(14)
 Polyomavirus(9)

Rash/Lesion


Monkeypox(2 - both clades)
Parvovirus(1)
 Rubella(1)
 Variola virus(1)

*** Note: Red Box indicates recent public health outbreaks world-wide**

Viral Surveillance Panel - Applications

U.S. declares monkeypox outbreak a public health emergency

By Helen Branswell Aug. 4, 2022



* **Sample Type:** Wastewater, Surfaces, Soil, Air, Clinical, Zoonotic, etc.

* We recommend ~1M reads for simple samples & ~5-10M reads for complex samples

Environmental Surveillance

- Emerging
- Outbreaks
- **Public Health**
- **USDA**
- **FDA**
- **BioDefense**
- **Customs Boarder Patrol**
- **Academic**
- **etc**

Clinical Surveillance

- Outbreaks
- Emerging
- Reflex
- **Public Health Labs**
- **Service Labs**
- **Government Labs**
- **Clinical Hospital Labs**
- **Academic Medical Centers**
- **etc**

Viral Metagenomics

- Research
- WGS
- **Academic**
- **Government**
- **etc**

NEW

Pan-Corona


For Research Use Only. Not for use in diagnostic procedures

illumina®

4

Introducing Pan-Coronavirus (Pan-CoV) Panel

- The Pan-CoV Panel allows for detection and **WGS of 203 known coronaviruses and over 370 strains of animal and closely related novel coronaviruses**, for surveillance of animal reservoir and sequencing of known human coronaviruses
- Hybrid capture panels **allow for up to 20% variation** in viral genome, allowing for mutations and novel strains to be detected
- At least one reference genome for each "branch" in the coronavirus phylogenetic tree is captured on the panel



Genera tiled end to end for genome enrichment and sequencing	Number of genus
Alphacoronavirus	69
Bafinivirus	1
Betacoronavirus	59
Deltacoronavirus	13
Gammacoronavirus	61
Grand Total	203

Illumina Pan-Coronavirus Panel (Pan-Cov Panel)

Genetically diverse coronavirus	Natural Host	Intermediate Host	Human Host
		?	→ HCoV-NL63
			→ HCoV-229E
			→ HCoV-OC43
		?	→ HCoV-HKU1
			→ SARS-CoV
			→ MERS-CoV
			→ SADS-CoV

Legend:
→ Spill over to intermediate hosts
→ Mild Infection
→ Severe Infection

Species	Number of genomes In Pan-Coronavirus panel from known hosts
Bat	64
Rodent	32
Mink, Pangolin, Murine, Porcine, snake	28
Bird	25
Bovine/Camel/Goat/Horse	11
Dolphin/whale/fish	9
Human	7
Dog	7
Feline/Ferret	6
Hedgehog	2

Pan-Coronavirus Panel - Applications

* **Sample Type:** Zoonotic, Wastewater, Surfaces, Air, etc.

* We recommend ~5-10M reads for complex samples



Environmental Surveillance

- Emerging
- Outbreaks
- **Zoonotic Researchers**
- **Public Health**
- **BioDefense**
- **Academic**
- **etc**

Viral Metagenomics

- Research
- WGS
- Discovery
- **Zoonotic Researchers**
- **Academic**
- **Government**
- **etc**



UPIP

For Research Use Only. Not for use in diagnostic procedures

illumina®

5

Urinary Pathogen Identification Panel (UPIP)

Powered by Explify™ UPIP

174

Pathogens



14
fungi



35
viruses



4
parasites



121
bacteria



46

Antibiotics

3,728

Antimicrobial Resistance Markers

Bonus: Simultaneous detection of STI pathogens

27 sexually transmitted pathogens detectable within the same urine specimen & workflow



Exact numbers of pathogens and AMR markers may change with product updates/improvements.

For Research Use Only. Not for use in diagnostic procedures.

Urinary Pathogen Identification Panel - Applications

* **Sample Type:** Urine, Swabs, Blood, Wastewater, etc.

* We recommend ~1M reads for simple samples & ~5-10M reads for complex samples

- UPIP provides **highly sensitive, culture-free identification and quantification** of both common and underrecognized uropathogens that can cause urinary tract infections.
- Also provides **Extensive antimicrobial resistance** information.



Clinical ID

- Chronic
- Undiagnosed
- Outbreaks
- **Public Health Labs**
- **Long Term Care Facilities**
- **Academic Medical Centers**
- **Commercial Service Labs**
- **Government Labs**
- **etc**

Clinical Surveillance

- Outbreaks
- Emerging
- Research
- Antimicrobial Resistance Monitoring (AMR)
- **Public Health Labs**
- **Government Labs**
- **Academic Research Centers**
- **etc**

Illumina Infectious Disease Panel Roadmap – summary comparison

	Illumina Respiratory Oligo Virus Panel RVOP	Illumina Respiratory Pathogen ID/AMR Panel RPIP	Illumina Urinary Pathogen ID/AMR Panel UPIP	Illumina Pan- Coronavirus Panel Pan-CoV	Illumina Viral Surveillance Panel VSP
Target(s)	40 common respiratory viruses	282 Respiratory pathogens (viral, bacteria, fungi) and 2108 antimicrobial resistance markers	174 Urinary pathogens (virus, bacteria, fungi, parasite) and 3728 antimicrobial resistance markers	225 Coronaviruses, over 370 strains	66 viral pathogens, 216 strains
Whole Genome vs. Detection	Whole Genome	Detection of pathogens and AMR markers Whole genome: SARS-CoV-2, Influenza A/B	Detection of pathogens and AMR markers	Whole Genome	Whole Genome
Library Prep	RNA Prep with Enrichment	RNA Prep with Enrichment	DNA Prep with Enrichment	RNA Prep with Enrichment	RNA Prep with Enrichment
Software	RNA Pathogen App	Explify RPIP	Explify UPIP	DRAGEN Microbial Enrichment App	DRAGEN Microbial Enrichment App
Configuration	Standalone panel	Panel, Library Prep, Indexes Set A, B , C, or D	Panel, Library Prep, Indexes Set A, B , C, or D	Standalone panel (96 reactions)	Standalone panel (96 reactions) Panel, Library Prep and Indexes Set A or B

Apps Galore!

1) Plug & Play, 2) Easy to Share & Collaborate, 3) Simple Push-Button Analysis, 4) Automation & Integration, 5) Secure

Illumina Core Apps



16S Metagenomics



GISAID
Submission



DRAGEN RNA
Pathogen Detection



DRAGEN COVID
Lineage

BaseSpace Labs Apps



Kraken
Metagenomics



Prokka Genome
Annotation



SRST2



Velvet de novo
Assembly



String Graph
Assembler



ITS
Metagenomics


BaseSpace 3rd Party Apps



DRAGEN
Metagenomics
Pipeline



MetaPhlAn



SPAdes Genome
Assembler



CosmosID
Metagenomics



EzBioCloud 16S-
based MTP



DNA Star



DeepChek-
HBV, HCV,
HIV




Explify RPIP
Analysis



SEAR



GENIUS
Metagenomics:
Know Now



QIIME

Breadth of ID and Surveillance @ illumina.....Future Panels?!

A world of use cases for Next-Generation Sequencing (NGS) for infectious disease



Antimicrobial resistance: Microbial sequencing, by NGS, is well established in research and public health settings and is emerging as a component that can be used to predict antimicrobial resistance. (Ranson E., 2020)



Early detection: NGS has been used in a vast array of analyses for infectious disease research of public health relevance to identify outbreak origins, track transmissions, investigate epidemic dynamics, determine etiological agents of a disease, and discover novel human pathogens (Berry, I. M., 2020)



Endemic Disease Surveillance: NGS is a great tool for vector control programs and public health environmental surveillance (Junqueira, A.C.M., 2017). That technology and bioinformatics tools offer new opportunities for Emerging Infectious Diseases surveillance in quality and in scale. (Kwok, K.T.T, 2020)




Wastewater: NGS can provide a comprehensive insight into microbial communities and their functional capacities in water and wastewater systems, thus eliminating the need to develop a new assay for each target organism or gene. (The Water Research Foundation, 2020)



Foodborne: Whole-Genome Sequencing provides information beyond the identity and relationship of strains; it can also help public health by improving the safety, quality, and shelf life of foods. (Brown, E., 2019)



Applying Genomics in Surveillance of Pathogens



Importance of Wastewater Surveillance for Public Health

Community-level monitoring of infectious diseases and detection of emerging pathogens



Asymptomatic population is included despite not seeking medical care/getting clinically tested



Track which variants are dominant in a community



Whole community-level view of infectious diseases regardless of access to healthcare



Discover emerging microbes in a population




Detect pathogens in community before hospitals can detect surges (pandemic preparedness)



Find antimicrobial resistance (AMR) genes

Wastewater Surveillance with Illumina NGS - Workflow

<https://www.illumina.com/destination/wastewater-surveillance-webinar-series.html>



3rd Party Solutions


- Health Care Acquired Infections
- TB - GenoScreen

Illumina Developed

Partnerships

Commercial Enablement

illumina®



Enablement through 3rd party solutions

Healthcare Acquired Infections (HAI)



BioMérieux

- Web based software solution to process bacterial whole genome sequence data (WGS).
- Detailed genetic characterization to accurately determine and track 13 of the most common bacteria found in healthcare settings.
- Uses illumina's WGS pipeline!

Mycobacterium tuberculosis (mtb)



GenoScreen

- An all-in-one assay based on targeted sequencing
- **DIRECT from Clinical Samples....no culture required!**
- IDs more than 140 mycobacterial species
- Prediction of resistance to 15 anti-TB drugs
- Uses illumina's Library Prep and Sequencers!

Enablement of our Solutions


- MiSeq
- NextSeq 1000/2000
- Infinity Long Reads
- New SBS Chemistry

Illumina Developed

Partnerships

Commercial Enablement

illumina®





MiSeq



Key Applications

- Bacterial & viral (small) genomes
- 16S profiling
- Targeted regions, amplicon panels



Foodborne Pathogen
Outbreaks



Public Health
England

Public Health Genomic
Epidemiology



Genomic Epidemiology



paired reads
Millions

kits available
Cycles

Nano

1

300

500

Micro

4

300

v2

15

50

300

500

v3

25

150

600



NextSeq 1000/2000 - will enable a broader range of applications and cost savings




Key Applications

- Larger genomes (i.e fungal)
- Metagenomics
- Whole Genome Sequencing at Scale
- Multiplexing many genomes
- Cost reduction via multiplexing

NEW


- **P1 and P2 with 600 cycle this year**
 - ~35hrs vs 65hrs (MiSeq)
 - ~same price, but 4x more data
- **P4 coming 2023/4 with new SBS Chemistry.**
 - Patterned flow cell like NovaSeq....future proof
- **NO Washes. Single, disposable cassette**
- **DRAGEN on-board for 6x faster secondary analysis**





Thank You!

Questions?



A Massive Leap in Sequencing Innovation with NovaSeq X

Powered by XLEAP-SBS™ Chemistry

Ultra high density
Flow cell

Advanced base
calling algorithms

Custom high-speed
camera sensor

Ambient shipped and
lyophilized reagents

Accelerated,
fully automated on
instrument analysis

Ultra high NA,
blue-green optics

