

LABORATORY STANDARD OPERATING PROCEDURE FOR WHOLE GENOME SEQUENCING ON MISEQ

Doc. No. PNL38

Ver. No. 03

Effective Date:

Page 1 of 31

1. **PURPOSE:** This SOP describes the standardized laboratory protocol for whole genome sequencing of bacterial organisms on the MiSeq, including instrument requirements and maintenance.
2. **SCOPE:** This document applies to all PulseNet WGS certified laboratories performing whole genome sequencing on enteric bacteria using the Illumina MiSeq platform and libraries generated with the library preparation kits accepted by PulseNet (Table 2) for submission of sequence data to PulseNet. Participating PulseNet laboratories may amend this procedure as necessary for use within their laboratories after validation per their laboratory's guidelines.
3. **DEFINITIONS**
 - 3.1. **BaseSpace:** Illumina cloud-based computing environment for next generation sequencing data analysis, management, and storage, including data sharing
 - 3.2. **CD:** Cluster Density; indicates the quantity of clusters that are generated per flow cell surface area during the cluster generation stage
 - 3.3. **CDC:** Centers for Disease Control and Prevention
 - 3.4. **Clusters PF:** clusters passing filter, percentage of generated clusters that pass an internal quality filtering procedure
 - 3.5. **CSV:** Comma-Separated Values (file) or Comma Delimited (file)
 - 3.6. **DNA:** Deoxyribonucleic acid
 - 3.7. **dsDNA:** double-stranded DNA
 - 3.8. **EBT:** Elution Buffer with Tween
 - 3.9. **Fastq:** A text-based file format for storing both a sequence and its corresponding quality scores
 - 3.10. **GB:** Giga byte
 - 3.11. **GHS:** Globally Harmonized System
 - 3.12. **HT1:** Hybridization buffer
 - 3.13. **IEM:** Illumina Experiment Manager
 - 3.14. **LRM:** Local Run Manager
 - 3.15. **MCS:** MiSeq Control Software
 - 3.16. **MSR:** MiSeq Reporter
 - 3.17. **NaOCl:** Sodium Hypochlorite
 - 3.18. **NaOH:** Sodium Hydroxide
 - 3.19. **nM:** Nanomolar
 - 3.20. **PhiX:** a control library that is derived from a well characterized bacteriophage genome. It has average size of 500 bp and balanced base composition.
 - 3.21. **PHL:** Public Health Laboratory
 - 3.22. **pM:** Picomolar
 - 3.23. **PN:** PulseNet
 - 3.24. **PPE:** Personal Protective Equipment
 - 3.25. **PulseNet Central:** PulseNet team at CDC comprising of the Outbreak Detection and Surveillance Unit and the Next Generation Subtyping Methods Unit
 - 3.26. **PulseNet/OutbreakNet SharePoint:** A closed, web-based collaboration application used for communication among PulseNet participants.
 - 3.27. **PR2:** Incorporation Buffer

LABORATORY STANDARD OPERATING PROCEDURE FOR WHOLE GENOME SEQUENCING ON MISEQ

Doc. No. PNL38

Ver. No. 03

Effective Date:

Page 2 of 31

- 3.28. **Q30:** The percentage of reads within the entire run that have Q scores >30.0
- 3.29. **QC:** Quality Control
- 3.30. **Q Score:** The sequence quality score for each individual base position in a sequence, indicating the accuracy of the base call. Phred scores are used, where $Q = -10\log(\text{Error Probability})$. The higher the quality score, the more reliable the base call. A Q30 means a 1 in 1000 likelihood of an incorrect base call at that position.
- 3.31. **RefID database:** a BioNumerics database, v 7.6 or higher, used for quality control of raw sequence data, assembly of sequences, contamination detection, and species identification. Part of the standard PulseNet workflow.
- 3.32. **RFID:** Radio-Frequency Identification
- 3.33. **RSB:** Resuspension Buffer
- 3.34. **SAV:** Sequencing Analysis Viewer, an application software that allows real-time viewing of quality metrics generated by the real-time analysis (RTA) software on the Illumina sequencing systems.
- 3.35. **SDS:** Safety Data Sheet
- 3.36. **SOP:** Standard Operating Procedure
- 3.37. **Tris-HCl:** Tris Hydrochloride
- 3.38. **UPS:** Uninterrupted Power Supply
- 3.39. **VP10:** Custom Read 1 Primer for sequencing Illumina DNA PCR-Free Prep libraries
- 3.40. **WGS:** Whole Genome Sequencing

4. RESPONSIBILITIES

4.1. PulseNet Public Health Laboratories:

- 4.1.1. Sequence isolates and perform quality check of the sequencing run and subsequent sequence data
- 4.1.2. Re-sequence any isolates that do not meet quality thresholds
- 4.1.3. Inform PulseNet Central, as necessary, about any complications with laboratory protocols or suspected instrument or reagent issues

4.2. PulseNet Central:

- 4.2.1. Perform additional sequence quality analysis in order to provide feedback and troubleshooting support for PHLs as necessary
- 4.2.2. Notify PN PHL if any sequences submitted do not meet quality thresholds
- 4.2.3. Communicate any suspected reagent issues to PHLs as necessary
- 4.2.4. Maintain and review SOPs on a regular basis and post on SharePoint

5. SAFETY

- 5.1. **Biosafety warning:** This document describes handling of DNA and associated products and does not describe best practices for handling of biological infectious material
- 5.2. **Chemical Safety Warning:** Take proper precautions and wear appropriate PPE when handling potentially hazardous chemicals. Ensure that chemicals, spent containers, and unused contents are disposed of in accordance with local and governmental safety standards. **See all relevant SDSs for additional information.**

LABORATORY STANDARD OPERATING PROCEDURE FOR WHOLE GENOME SEQUENCING ON MISEQ

Doc. No. PNL38

Ver. No. 03

Effective Date:

Page 3 of 31

- 5.2.1. The MiSeq reagent cartridges contain formamide (GHS classification Category 1B for reproductive toxicity), an aliphatic amide that is a potential reproductive toxin. Personal injury can occur through inhalation, ingestion, skin and eye contact.
- 5.2.2. Ethanol is flammable (GHS Flammability Category 2)
- 5.2.3. Sodium hydroxide is corrosive (GHS Category 1A and 1, GHS Category 3 for acute hazards to aquatic environment)

6. REAGENTS

6.1. MiSeq Reagent Kit Options:

- v2 Nano, 300 Cycles – Illumina, Cat# MS-103-1001
- v2 Nano, 500 Cycles – Illumina, Cat# MS-103-1003
- v2 Micro, 300 Cycles – Illumina, Cat# MS-103-1002
- v2 300 Cycles – Illumina, Cat# MS-102-2002
- v2 500 Cycles – Illumina, Cat# MS-102-2003
- v3 600 Cycles – Illumina, Cat# MS-102-3003

NOTE: *PulseNet does not accept in any circumstances sequencing chemistries below 300 cycles.*

6.1.1. Box 1 of 2. Store at -15°C to -25°C, away from light.

- MiSeq Cartridge (see Section 5.2.1 for chemical safety information)
- Hybridization Buffer (HT1)

6.1.2. Box 2 of 2. Store at 2-8°C.

- Incorporation Buffer (PR2)
- Flow Cell

6.2. Ethanol, molecular-grade, 95-100% (Fisher Scientific, Cat# BP2818-500 or equivalent)

6.3. Ethanol, lab-grade, 70% or equivalent for disinfection purposes (Fisher Scientific, Cat# 04-355-309 or equivalent)

6.4. Water, molecular grade (Fisher Scientific, Cat# BP24701 or equivalent)

6.5. Sodium Hypochlorite, lab-grade, 4 – 8.25% (Fisher Scientific Cat# SS290-1 or equivalent).
NOTE: *It is recommended to use laboratory grade sodium hypochlorite, due to possible additives in over-the-counter grade sodium hypochlorite, which may interfere with the instrument's fluidics*

6.6. Sodium Hydroxide, suitable for cell culture, 1N (Millipore Sigma, Cat# S2770-100ML or equivalent), pH should be ≥ 12.5 .

NOTE: *Recommended to use liquid over powdered. Also advisable to aliquot and freeze in 100 μ l, single-use aliquots.*

6.7. Tween 20, molecular grade, (Millipore Sigma, Cat# P9416-50ML or equivalent)

6.8. **OPTIONAL:** 10nM PhiX Control Kit v3 (Illumina, Cat# FC-110-3001). Store at -25°C to -15°C

6.9. **If sequencing Illumina DNA PCR-Free Prep libraries:** Illumina DNA PCR-Free Prep Sequencing Primers Read 1 (Illumina, Cat. No. 20041796, 4 x 7.5 ml)

7. SUPPLIES

7.1. Ice

7.2. Illumina MiSeq Wash Tray (n=2, one for bleach washes, another for water washes only)

7.3. Lint-free wipes (Fisher Scientific, Cat# 06-666 or equivalent)

LABORATORY STANDARD OPERATING PROCEDURE FOR WHOLE GENOME SEQUENCING ON MISEQ

Doc. No. PNL38

Ver. No. 03

Effective Date:

Page 4 of 31

- 7.4. Lens paper (Fisher Scientific, Cat# 11-996 or equivalent)
- 7.5. MiSeq Disposable Wash Tube (Illumina, Cat# MS-102-9999)
- 7.6. Wash Buffer Bottle for MiSeq, Nalgene container, 500 ml (n=2, for 0.5% Tween 20 and molecular grade water, Fisher Scientific, Cat# 15-350-205 or equivalent).
ALTERNATIVELY: Can use spent PR2 bottles.
- 7.7. Microseal B film (BioRad, Cat# MSB-1001)
- 7.8. 1.5 ml microcentrifuge tubes (Fisher Scientific, Cat# 05-408-129 or equivalent)
- 7.9. Sterile pipette tips, filtered: 20 µl, 200 µl, 1000 µl and 1000 µl extended length volumes (Rainin, Cat# 30389225, 30389239, 30389212 and 30389223 or equivalent)
- 7.10. Wash bottle with nozzle (x3, one each for 0.5% Tween 20, 100% ethanol, and molecular grade water (Fisher Scientific, Cat# FB0340922D or equivalent))

8. EQUIPMENT

- 8.1. MiSeq Desktop Sequencer
- 8.2. Micropipettes, capable of volumes from 1 µl to 1000 µl
- 8.3. Heat block
- 8.4. Ice buckets
- 8.5. Pipet-aid
- 8.6. Microcentrifuge
- 8.7. Microplate centrifuge or equivalent
- 8.8. **OPTIONAL:** UPS back-up for the MiSeq (Recommended: Staples, Cyberpower AVR Series Line Interactive 1.5 kVA UPS, Cat# CP1500AVRLCD)
- 8.9. **OPTIONAL:** External encrypted hard drive or server for data transfer and storage if BaseSpace or networking of the instrument is not available (Suggested: CDW, DataLocker H350 Basic Hard Drive 1 TB USB 3.0, Cat# 4075102)

9. PROCEDURE

NOTE: It has been documented and observed that cluster density can be negatively affected by the use of ammonium-based cleaning products near sequencing equipment, including lab benches and pipets used for library preparation. Do not use quaternary ammonium compounds or wipes near or on sequencing equipment!

9.1. Prepare the Reagent Cartridge for the MiSeq

- 9.1.1. Plan sufficient time for thawing the sequencing reagents (sequencing cartridge, HT1 buffer) and the V10 primer (if sequencing Illumina DNA PCR-Free Prep libraries). The HT1 and the VP10 primer buffer may be thawed on ice or at 2-8°C. See Table 1 for guidance regarding thaw time, method and storage for the sequencing cartridge:

Method of Thawing	Time Required to Thaw	Thawed Cartridge Storage Life (2 – 8°C or on ice)	A cartridge should not be re-frozen after thawing
Room temperature water bath (opened)	~ 1 hour	24 hours ¹	
2 – 8°C (unopened)	8-12 hours	7 days	

LABORATORY STANDARD OPERATING PROCEDURE FOR WHOLE GENOME SEQUENCING ON MISEQ

Doc. No. PNL38

Ver. No. 03

Effective Date:

Page 5 of 31

Table 1. Guidance for reagent kit thawing and storage after thawing.

¹Illumina recommendation 6 h

- 9.1.2. Once thawed, invert the reagent cartridge 10 times to mix, and then visually inspect to verify that all positions are completely thawed and are free of precipitates
- 9.1.3. Gently tap the cartridge on the bench to reduce air bubbles in the bottom of the reagent wells and ensure that the foil seals covering the wells are unobstructed
- 9.1.4. Place the thawed reagent cartridge and HT1 on ice or set aside at 2 – 8°C until the sample is ready to be loaded

9.2. Chemical Denaturation of the Pooled Library DNA

NOTE: Refer to the available PulseNet Library Preparation SOPs for various kit options in prepping and pooling the libraries. If using the DNA Prep Checklist (PNL35.W3), pooling, denaturation and dilution and instrument loading are listed as steps #73-97.

- 9.2.1. Pre-heat a heat block to 96°C ± 1.

NOTE: Skip the heat denaturation step if sequencing the Illumina DNA PCR-Free Prep libraries.

- 9.2.2. Retrieve the pooled library plate/tubes, and if frozen, thaw on ice. Centrifuge at 800 – 1200 rpm (or 280 x g) for 1 minute.
- 9.2.3. Remove a 100 µl aliquot of 1 N NaOH from the freezer, and thaw on ice
- 9.2.4. Dilute to 0.2 N by adding 400 µl of molecular grade water to the 100 µl aliquot of 1 N NaOH. Invert the tube several times to mix.

NOTE: Freshly diluted NaOH should be used within 6 hours of dilution.

- 9.2.5. Combine 5 µl of the pooled sample library and 5 µl of 0.2 N NaOH in a new 1.5 ml microcentrifuge tube (now at 1 or 2 nM)
- 9.2.6. Incubate for 5 minutes at room temperature to denature the dsDNA
- 9.2.7. **IMMEDIATELY** add 990 µl of pre-chilled HT1 to the tube containing 10 µl of denatured pooled library after the incubation. Refer to Table 2 for initial and denatured library concentrations.

Library Prep Kit	Initial Library Pool Concentration	Denatured Library Concentration After Addition of 990 µl of HT1
Nextera XT, NEB Ultra II FS, QIAseq FX, & KAPA HyperPlus	2 nM	10 pM
Illumina DNA Prep, Nextera XT, NEB Ultra II FS, QIAseq FX, & KAPA HyperPlus	4 nM	20 pM
Illumina DNA PCR-Free	2 nM	20 pM

Table 2. Denatured library concentration from initial library pool concentration

9.3. Dilution to Desired Final Loading Concentration and Heat Denaturation of the Pooled Library DNA

LABORATORY STANDARD OPERATING PROCEDURE FOR WHOLE GENOME SEQUENCING ON MISEQ

Doc. No. PNL38

Ver. No. 03

Effective Date:

Page 6 of 31

- 9.3.1. Determine desired final loading concentration, based on library preparation type and sequencing kit to be used. See Table 3 below for guidance.

NOTE: Final loading concentration may also be affected by a particular instrument and/or user and may be adjusted above or below the recommended concentrations in Table 3 in order to achieve optimal cluster density (CD).

Library Prep Kit	Recommended Final Loading Concentration	
	v2	v3
Nextera XT	10 pM	15 pM
Illumina DNA Prep (DNA Flex)	12 pM	15 pM
NEB Ultra II FS & KAPA HyperPlus	10 pM	ND ¹
QIAseq FX	8 pM	ND ¹
Illumina DNA PCR-Free	12 pM	ND ¹

Table 3. Recommended loading concentrations for different library prep kits. ¹Not determined

- 9.3.2. Refer to Table 4 below to obtain the final loading concentration of the pooled library and dilute accordingly:

Final Loading Concentration	Denatured Pool Concentration = 10 pM		Denatured Pool Concentration = 20 pM	
	Required Volume of HT1	Required Volume of Denatured Pool	Required Volume of HT1	Required Volume of Denatured Pool
8 pM	200 µl	800 µl	600 µl	400 µl
10 pM	NA	NA	500 µl	500 µl
12 pM	NA	NA	400 µl	600 µl
15 pM	NA	NA	250 µl	750 µl

Table 4. Required pool and HT1 volumes for dilution to desired final concentration

- 9.3.3. Thoroughly mix the diluted libraries.

OPTIONAL: *PhiX control may be spiked in at this point. See Appendix PNL38-1 for instructions and more information.*

- 9.3.4. Heat the denatured, diluted library pool to $96 \pm 1^\circ\text{C}$ for 2 minutes in a heat block to ensure complete denaturation.

LABORATORY STANDARD OPERATING PROCEDURE FOR WHOLE GENOME SEQUENCING ON MISEQ			
Doc. No. PNL38	Ver. No. 03	Effective Date:	Page 7 of 31

NOTE: Skip the heat denaturation step (9.3.4-9.3.5) if sequencing the Illumina DNA PCR-Free Prep libraries.

- 9.3.5. Immediately cool in an ice-water bath (made by combining approximately 3 parts ice and 1 part water) for at least 5 minutes prior to loading.

NOTE: The DNA library may sit on ice or at 2 – 8°C until ready for loading (< 30 minutes). If greater than 30 minutes passes before starting the run, it may be necessary to re-denature and dilute the pool.

9.4. **Prepare the MiSeq Instrument for the Run**

- 9.4.1. Prior to setting up the run, ensure that there is enough free disk space (100 GB) on the instrument. If there is less than 100 GB available, see Section 9.11.3. for instructions on deleting files prior to loading the instrument and starting a run.

- 9.4.2. Ensure that the Sample Sheet (and Sample Plate, if applicable) setup has been completed in the Local Run Manager (LRM) or Illumina Experiment Manager (IEM) on the instrument. See Appendixes PNL38-2 and PNL38-3, respectively, for instructions.

NOTE: The workbooks contain separate tabs for different software versions on the MiSeq and additional Illumina instruments: IEM, MiSeq LRM3, MiSeq LRM versions earlier than 3, iSeq LRM, and MiniSeq LRM. When exporting the sample sheet, ensure that the proper tab is used for the software available on the MiSeq.

- 9.4.3. Open the MiSeq Control Software window, select “Sequence”

- 9.4.3.1. **If using IEM**, proceed immediately to Preparing and Loading the Flow Cell.

NOTE: If using BaseSpace, on the following screen, check the box to make sure the user info fields are active and enter credentials to ensure the run data will be streamed to BaseSpace.

- 9.4.3.2. **If using LRM**, select “Local Run Manager”, select appropriate BaseSpace streaming option, and confirm run parameters before selecting “Next”

9.5. **Prepare and Load the Flow Cell**

- 9.5.1. Retrieve the flow cell and the Incorporation buffer (PR2) from storage (2 – 8°C).

NOTE: Do not keep or store the PR2 buffer on ice as it may form particulates.

- 9.5.2. While wearing clean, powder-free gloves, carefully remove the flow cell from the container, without touching the glass of the flow cell

- 9.5.3. Rinse the flow cell assembly with molecular-grade water, making sure that both the glass and plastic casing are thoroughly rinsed of excess salts

- 9.5.4. Pat the flow cell dry with a lint-free wipe, using care around the black port gasket. Make sure to remove all excess fluid.

- 9.5.5. Wet a clean piece of lens paper with ethanol and clean the flow cell glass, making sure that the glass imaging area (B in Fig. 1 below) is free of streaks, lint and tissue fibers.

NOTE: Do not add ethanol directly to flow cell. Avoid getting ethanol on the flow cell gaskets (A & C in Fig. 1).

LABORATORY STANDARD OPERATING PROCEDURE FOR WHOLE GENOME SEQUENCING ON MISEQ

Doc. No. PNL38

Ver. No. 03

Effective Date:

Page 8 of 31




Figure 1. Flow cell. Port gaskets (A & C), imaging channels (B)

- 9.5.6. Dry any excess alcohol with the lens paper, and visually inspect to make sure that the flow cell ports (A & C) are free of obstructions, that the gasket is well-seated around the flow cell ports, and that there are no smudges or debris on the glass
- 9.5.7. Raise the flow cell compartment door on the MiSeq and press the release button to open the clamp and remove the previously used flow cell
- 9.5.8. Ensure that the flow cell stage is free of lint. If necessary, wipe with an alcohol wipe or similar and allow to dry.
- 9.5.9. Place the flow cell on the stage, close the flow cell clamp (it will click when closed fully), close the compartment door and select “Next” in the MiSeq Control Software

9.6. Load the Incorporation Buffer

- 9.6.1. Gently invert the Incorporation buffer (PR2) to mix, and then remove the lid
- 9.6.2. Open the reagent compartment door and raise the sipper handle until it locks into place
- 9.6.3. Remove the Wash Buffer bottle, re-cap, and set aside for future use
- 9.6.4. Place the PR2 bottle where the MiSeq Wash Buffer bottle was previously located
- 9.6.5. Remove the waste bottle and empty into an appropriate container for chemical waste, if necessary.

NOTE: All liquid that accumulates in the MiSeq waste bottle should be disposed of as formamide waste, including liquid from wash cycles. Ensure that chemical waste is labeled and disposed of appropriately.

- 9.6.6. Slowly lower the sipper handle. Make sure that the sippers lower completely into the PR2 and waste bottles.
- 9.6.7. Check the lower-right corner of the screen to confirm that the RFID of the PR2 bottle was read successfully.

NOTE: If the RFID is not read by the system for any step, the software prompts you through the steps to obtain a temporary bypass code and proceed with setting up the run. For more information, see “Resolving RFID Read Failure” in the MiSeq System User Guide.

- 9.6.8. Select “Next” on MCS

9.7. Prepare and Load the Cartridge

- 9.7.1. Remove the thawed reagent cartridge from storage, dry it thoroughly, and mix the reagents by inverting gently 10 times.

LABORATORY STANDARD OPERATING PROCEDURE FOR WHOLE GENOME SEQUENCING ON MISEQ

Doc. No. PNL38

Ver. No. 03

Effective Date:

Page 9 of 31

NOTE: *Mixing of reagents is very important! If the library pool is put into the cartridge prior to mixing the cartridge reagents, use an extended pipet tip to remove it, tape the sample reservoir, mix the cartridge and then re-load the pool.*

- 9.7.2. Tap the reagent cartridge on a hard surface to collect contents at the bottom of the reservoirs
- 9.7.3. Use a clean 1000 µl pipette tip to pierce the foil seal over the reservoir labeled “Load Sample”
- 9.7.4. Load **600 µl** of the denatured DNA library into the “**Load Sample**” reservoir (Fig. 2).
NOTE: *If sequencing Illumina DNA PCR-Free Prep libraries, add 600 µl of the thawed VP10 primer into the custom reservoir, position 18 (A in Fig. 2).*




Figure 2. Reserved reservoirs in the sequencing cartridge: No. 17 for the library pool to be sequenced, no. 18-20 for custom primers.

NOTE: *It is recommended to load the sample and primer toward the bottom of reservoir using an extended 1000 µl pipette tip to reduce the risk of droplets forming on the sides.*

Take care to avoid droplets on the foil seal as the sample is dispensed.

- 9.7.5. Tap gently on a hard surface to bring all contents to the bottom
- 9.7.6. Open the reagent chiller door and remove the wash tray. If necessary, dry the bottom of the chiller compartment using absorbent wipes.
NOTE1: *Do NOT forcefully pull the wash tray; if there is resistance, the sippers may still be lowered. Wait until the MiSeq gives the prompt to load the cartridge and try again.*
NOTE2: *Do not leave the reagent chiller door open for extended periods of time.*
- 9.7.7. Slide the reagent cartridge into the chiller compartment until it stops
- 9.7.8. Close the chiller door and check the screen to confirm that the RFID of the reagent cartridge has been read successfully
- 9.7.9. Close the reagent compartment door. **If using LRM**, proceed to step 9.7.11.
- 9.7.10. **If using IEM:**
 - 9.7.10.1. Select “Change Sample Sheet” and in the subsequent screen, select the “Browse” button, and navigate to the appropriate sample sheet
 - 9.7.10.2. Select the desired sample sheet and click “OK”
 - 9.7.10.3. Select “Next” to proceed to the run parameters review screen
- 9.7.11. Review Run Parameters and ensure that they are correct (Experiment Name, Analysis Workflow, Read Length, etc.) then select “Next” to proceed to the Pre-Run Check. **If**

LABORATORY STANDARD OPERATING PROCEDURE FOR WHOLE GENOME SEQUENCING ON MISEQ

Doc. No. PNL38

Ver. No. 03

Effective Date:

Page 10 of 31

using LRM, proceed to step 9.8 as the run will start automatically once pre-run checks complete successfully.

NOTE: The system performs a check of all run components, disk space, and network connections prior to run start. If any part of the pre-run check fails, a message will appear on the screen with general instructions describing the error or detailing how to correct it. For more information, see “Resolving Run Setup Errors” in the MiSeq System User Guide. If all items successfully pass the pre-run check, the system is ready to start the run.

9.7.12. Select “Start Run”.

NOTE1: Once the run has been started, do not open the flow cell compartment or the reagent compartment doors, neither should the instrument monitor be engaged unless the run is to be stopped or paused.

NOTE2: Image capture on the MiSeq is sensitive to vibration. Performing tasks that cause vibration near/on the instrument during a run could cause the run to halt or adversely impact sequencing results.

9.8. Review Run Metrics

9.8.1. Upon run completion, confirm that the sequencing run meets basic quality metrics, which will be shown on MCS when the run is complete (Fig. 3)




Figure 3. MiSeq Control Software Run Completion Screen. Primary metrics are circled in red.

NOTE: These run metrics may also be found on the Summary Tab of SAV, see Fig. 4 below.

LABORATORY STANDARD OPERATING PROCEDURE FOR WHOLE GENOME SEQUENCING ON MISEQ

Doc. No. PNL38

Ver. No. 03

Effective Date:

Page 11 of 31

Analysis | Imaging | Summary | Indexing

Cycle Status

Extracted: 518 Called: 518 Scored: 518 Error Rated: 0

Run Summary

Level	Yield Total (G)	Projected Total Yield (G)	Aligned (%)	Error Rate (%)	Intensity Cycle 1	% >= Q30
Read 1	3.62	3.62	0.00	NaN	172	94.31
Read 2 (I)	0.10	0.10	0.00	NaN	851	97.90
Read 3 (I)	0.10	0.10	0.00	NaN	445	96.10
Read 4	3.62	3.62	0.00	NaN	126	82.74
Non-indexed Total	7.45	7.24	0.00	NaN	149	88.59
Total	7.45	7.45	0.00	NaN	399	88.76

Read 1

Lane	Tiles	Density (K/mm ²)	Cluster PF (%)	Legacy Phasing/Prephasing Rate	Phasing slope/offset	Prephasing slope/offset	Cluster Count Raw (M)	Cluster Count PF (M)	% >= Q30	Yield (G)	Cycles Err Rated	Aligned (%)	Error Rate (%)	Error Rate 35 cycle (%)	Error Rate 75 cycle (%)
1	28	781 ± 18	95.06 ± 0.26	0.059 / 0.048	NaN / NaN	NaN / NaN	15.24	14.49	94.31	3.62	0	0.00 ± 0.00	NaN ± NaN	NaN ± NaN	NaN ± NaN

Read 2 (I)

Lane	Tiles	Density (K/mm ²)	Cluster PF (%)	Legacy Phasing/Prephasing Rate	Phasing slope/offset	Prephasing slope/offset	Cluster Count Raw (M)	Cluster Count PF (M)	% >= Q30	Yield (G)	Cycles Err Rated	Aligned (%)	Error Rate (%)	Error Rate 35 cycle (%)	Error Rate 75 cycle (%)
1	28	781 ± 18	95.06 ± 0.26	0.000 / 0.000	NaN / NaN	NaN / NaN	15.24	14.49	97.90	0.10	0	NaN ± NaN	NaN ± NaN	NaN ± NaN	NaN ± NaN

Read 3 (I)

Lane	Tiles	Density (K/mm ²)	Cluster PF (%)	Legacy Phasing/Prephasing Rate	Phasing slope/offset	Prephasing slope/offset	Cluster Count Raw (M)	Cluster Count PF (M)	% >= Q30	Yield (G)	Cycles Err Rated	Aligned (%)	Error Rate (%)	Error Rate 35 cycle (%)	Error Rate 75 cycle (%)
1	28	781 ± 18	95.06 ± 0.26	0.000 / 0.000	NaN / NaN	NaN / NaN	15.24	14.49	96.10	0.10	0	NaN ± NaN	NaN ± NaN	NaN ± NaN	NaN ± NaN

Read 4

Lane	Tiles	Density (K/mm ²)	Cluster PF (%)	Legacy Phasing/Prephasing Rate	Phasing slope/offset	Prephasing slope/offset	Cluster Count Raw (M)	Cluster Count PF (M)	% >= Q30	Yield (G)	Cycles Err Rated	Aligned (%)	Error Rate (%)	Error Rate 35 cycle (%)	Error Rate 75 cycle (%)
1	28	781 ± 18	95.06 ± 0.26	0.103 / 0.043	NaN / NaN	NaN / NaN	15.24	14.49	82.74	3.62	0	0.00 ± 0.00	NaN ± NaN	NaN ± NaN	NaN ± NaN

Copy to Clipboard

Zip My Run

Figure 4. Screenshot of the SAV Summary Tab highlighting the key run metrics to be reviewed

- 9.8.2. Record the primary run metrics on the workbook. Compare Q30, cluster density and clusters passing filter to PulseNet's recommended run metrics thresholds in Table 5.

NOTE: If the run metrics do not meet the thresholds listed in Table 5, it is possible that sequences may pass the critical quality metrics outlined in PNQ07. Further analysis of the sequences in the RefID database (PND20) is required to determine whether individual sequences are acceptable or need to be repeated.

Kit Chemistry	Q30 (%)	Cluster Density (K/mm ²)	Clusters Passing Filter (%)
v3, 500 cycle*	≥ 70	1200-1400	~ 80 or higher
v2, 500 cycle	≥ 75	600-1200	
v2, 300 cycle	≥ 80	600-1200	
Nano, v2 300 cycle	≥ 80	600-1200	
Nano, v2 500 cycle	≥ 75	600-1200	
Micro, v2 300 cycle	≥ 80	600-1200	

Table 5. Run metric thresholds for sequencing chemistries accepted by PulseNet. *Kit standard is 600c but should be sequenced at 500 cycles for better quality as suggested by PulseNet.

9.9. Export Data from the MiSeq

NOTE1: This section is only necessary for users unable to stream to BaseSpace or if a secondary data repository (e.g., local server) is desired for local data analysis or storage.

LABORATORY STANDARD OPERATING PROCEDURE FOR WHOLE GENOME SEQUENCING ON MISEQ

Doc. No. PNL38

Ver. No. 03

Effective Date:

Page 12 of 31

NOTE2: There will be 2 fastq files per isolate for paired reads (R1 and R2). This data may be copied onto an external hard drive and moved to a computer workstation. Alternatively, the MiSeq instrument can also be connected to a local network share drive to directly transfer files from the instrument.

9.9.1. **For IEM,** data is found on the D:\ drive → Illumina → MiSeqOutput → Run folder (choose the most recent run folder) → Data → Intensities → BaseCalls. Within the BaseCalls folder, the raw read data is in fastq.gz format.

9.9.2. **For LRM,** data is found on the D:\ drive → Run folder (choose the most recent run folder) → Alignment_1 → Run folder → Fastq.gz)

NOTE1: It may take up to 1 hour for the MiSeq to generate fastq files post-run as the Real-Time Analysis Software lags behind the sequencing run. If using BaseSpace, the data can be accessed through the run or in a designated project.

NOTE2: If fastq data is unavailable for a given isolate(s) due to incorrect index assignment or sample sheet error, the run can be requeued for analysis using either MiSeq Reporter (MSR), LRM or corrective options in BaseSpace (see Appendix PNL38-4).

9.10. Post-Run Washes and Reagent Disposal

9.10.1. Once the run is complete, select the “Start Wash” option at the bottom right-hand side of the screen

9.10.2. Ensure the “Perform optional template line wash” checkbox is checked to proceed with a Post-Run bleach wash. Sodium hypochlorite is used to reduce carry-over of nucleic acid from previous runs.

9.10.3. Prepare a fresh dilution of molecular-grade sodium hypochlorite (NaOCl) using molecular-grade water. Refer to Table 6 below for serial dilution steps. Choose the first dilution according to the concentration of stock NaOCl.

NOTE1: The second dilution (0.01%) may be prepared directly in a MiSeq wash tube/reservoir (A in Fig. 2).

NOTE2: Use of regular household bleach has not been validated and is discouraged.

Dilution	Starting NaOCl Concentration	NaOCl Volume	Water Volume	Final NaOCl Concentration
First Dilution:	4%	30 µl	570 µl	0.20%
	5%	24 µl	576 µl	
	6%	20 µl	580 µl	
	8.25%	15 µl	585 µl	
Second Dilution:	0.20%	50 µl	950 µl	0.01%

Table 6. The calculations needed to prepare 0.2% and 0.01% dilutions from stock NaOCl

9.10.4. Insert the wash tube into position 17 in the designated bleach wash tray until the neck is flush with the tray (B in Fig. 5).

LABORATORY STANDARD OPERATING PROCEDURE FOR WHOLE GENOME SEQUENCING ON MISEQ

Doc. No. PNL38

Ver. No. 03

Effective Date:

Page 13 of 31

NOTE: It is recommended to have 2 wash trays, one designated for bleach and Tween 20 and the other designated for water wash only.




Figure 5. Inserting the 0.01% Sodium hypochlorite wash tube in wash tray

9.10.5. Fill the remaining reservoirs with 0.5% Tween 20 wash solution.

NOTE: 0.5% Tween 20 can be prepared by adding 5 ml of Tween 20 into 995 ml of molecular-grade water. Only prepare 1 L of 0.5% Tween 20 at a time to ensure its use in a timely fashion. Over time, the detergent may separate out of solution.

9.10.6. Remove excess moisture from the surface of the wash cartridge with a lint-free wipe

9.10.7. Open the reagent compartment door and chiller door and remove the used cartridge. Set the cartridge aside for disposal.

NOTE: Reagent well 8 of the cartridge contains formamide and must be disposed of as hazardous chemical waste. See Section 5 and the SDS for safety details.

9.10.8. Slide the wash tray into the chiller compartment until it stops, then close the chiller door

9.10.9. Raise the sipper handle, remove the PR2 bottle and replace with a MiSeq Wash Buffer bottle containing 300-500 ml of 0.5% Tween 20 solution.

NOTE: Replace the Wash Buffer bottle as needed to limit contaminants.

9.10.10. Remove the waste container and dispose of contents following the laboratory's hazardous chemical waste disposal methods (See Section 5.2.1. for safety information).

NOTE: All liquid that accumulates in the MiSeq waste bottle should be disposed of as toxic formamide waste, including liquid from wash cycles. All formamide waste bottles must be clearly labeled as chemical waste and disposed of appropriately.

9.10.11. Return the emptied waste container to the reagent compartment, lower the sipper arm, and close the reagent compartment door

9.10.12. Select "Next" to begin the first post-run template line wash. This will take approximately 30 minutes.

9.10.13. Once the post-run template line wash is complete, select "Done" on the screen to return to the MiSeq Control Software home screen

9.10.14. Immediately following the template line wash, a second Post-Run wash with molecular-grade water must be completed to ensure any trace amount of sodium hypochlorite is removed from the template line. From the home screen select "Perform Wash" > "Perform Post-Run Wash".

LABORATORY STANDARD OPERATING PROCEDURE FOR WHOLE GENOME SEQUENCING ON MISEQ

Doc. No. PNL38

Ver. No. 03

Effective Date:

Page 14 of 31

9.10.15. Select “Start Wash”.

NOTE: The “Template Line Wash” box should remain **un-checked**.

9.10.16. Fill all reservoirs of the designated water wash tray with molecular-grade water

9.10.17. Remove excess moisture from the surface of the wash cartridge with a lint-free wipe

9.10.18. Remove the wash tray containing Tween 20 and the bleach wash tube from the chiller compartment and rinse the wash tube thoroughly. Rinse the wash tray with de-ionized water (or equivalent) and let dry.

9.10.19. Place the water wash tray in the chiller compartment and close the door

9.10.20. Raise the reagent sipper arm and replace the 0.5% Tween 20 wash bottle with a wash bottle containing 300-500 ml of molecular-grade water and lower the sipper arm

9.10.21. Close the reagent compartment door and select “Next” to begin the wash

9.10.22. Once the wash is complete, select “Done” to complete the wash cycle and return to the home screen. Keep the flow cell, wash bottle, wash cartridge and waste receptacle in place until the next sequencing run is to be performed

9.11. Instrument Maintenance

	Maintenance Frequency		
	Following A Run	Weekly	Monthly
Type of Maintenance	Post-Run Template Line Wash, followed by Post-Run Wash (using water only)	Post-Run Wash and Reboot	Maintenance Wash + Power cycle
Time Required	30 minutes + 20 minutes	20 minutes + 10 minutes	60 minutes+ 20 minutes
Reagents Needed	- 0.5% Tween 20 - 0.01% NaOCl - Molecular grade water	- 0.5% Tween 20	- 0.5% Tween 20

Table 7. Summary of the maintenance of the MiSeq

9.11.1. Weekly Maintenance:

9.11.1.1. Wash: A Post-Run Wash with 0.5% Tween 20 should be performed to ensure proper instrument fluidics performance. If a sequencing run has been performed that week, the routine post-run wash followed by the water wash will suffice as the weekly wash (Step 9.10.). If the instrument has not been used, a manual post-run wash using 0.5% Tween 20 should be performed.

9.11.1.2. Reboot: The instrument should be rebooted weekly and after data deletion when freeing the virtual memory.

NOTE1: The weekly reboot can be skipped on the week when the monthly power cycle (step 9.11.2.2.) has been performed.

NOTE2: When performing the reboot, having an external hard drive plugged in may cause an error at the start-up.

LABORATORY STANDARD OPERATING PROCEDURE FOR WHOLE GENOME SEQUENCING ON MISEQ

Doc. No. PNL38

Ver. No. 03

Effective Date:

Page 15 of 31

- 9.11.1.2.1. On the MCS home screen, select “Manage Instrument” then “Reboot”. It will take approximately 10 minutes for the system to reboot and start the MiSeq Control software.
- 9.11.2. Monthly Maintenance:
 - 9.11.2.1. Maintenance Wash: To be performed at least once a month. This is a series of 3 washes that may require refilling of the wash tray after each wash cycle (20 minutes each). The total time for a maintenance wash is approximately 60 minutes.
 - 9.11.2.1.1. Under Wash Options, select “Maintenance Wash”
 - 9.11.2.1.2. Fill the wash tray and wash buffer bottle with 0.5% Tween 20
 - 9.11.2.1.3. After the first two washes, if prompted, top off the reservoirs of the wash tray and wash buffer bottle with 0.5% Tween 20.
NOTE: *If desired, the third wash can be a water wash.*
 - 9.11.2.2. Power cycle:
 - 9.11.2.2.1. Select the “Manage Instrument” menu option on the home screen and then the “Shut Down” option
 - 9.11.2.2.2. Once the instrument has shut down (noted by a distinct click sound), reach behind the right-hand side of the instrument and flip off the power switch located near the power cord
 - 9.11.2.2.3. Let the instrument power down for at least 5 minutes and then flip the switch back on. The instrument will boot up the MCS in about 10 minutes.
- 9.11.3. Management of Disk Space on the Data (D:\) Drive:
 - 9.11.3.1. The MiSeq requires 100 GB of free disk space on the D:\ drive before starting a sequencing run. **Data should be cleared only *after* any data transfer/back-up is complete, and if no troubleshooting is necessary for that run.** The following instructions outline which files may be deleted and their locations:
 - 9.11.3.1.1. Navigate to Data (D:\ Drive)/Illumina/MiSeqOutput/(most recent run folder)/. Within this folder select the “Images” and/or “Thumbnail_images” folders and delete; suggested for runs that are over one month old.
 - 9.11.3.1.2. Navigate to Data (D:\ Drive)/Illumina/MiSeqAnalysis/. Delete any run folders excluding the most recent folder.
 - 9.11.3.1.3. Empty the Recycle Bin
 - 9.11.3.1.4. Reboot the system following large data deletions
 - 9.11.3.2. **BaseSpace Users:** It is recommended to maintain a duplicate copy of the run folder locally on the instrument as a backup. Doing the following will ensure that a run folder will be generated in the Data drive file folders, in addition to being streamed to BaseSpace:
 - 9.11.3.2.1. From the MCS home screen, select “Run Options”
 - 9.11.3.2.2. Under the “Run Settings” tab, check the “When using BaseSpace, replicate analysis locally on MiSeq” checkbox
 - 9.11.3.2.3. Select “Save and Return” to save these changes and return to the home screen

LABORATORY STANDARD OPERATING PROCEDURE FOR WHOLE GENOME SEQUENCING ON MISEQ

Doc. No. PNL38

Ver. No. 03


Effective Date:

Page 16 of 31


NOTE: *This run folder may be deleted if run data was successfully streamed to BaseSpace and local data storage is not required.*

10. FLOW CHARTS:

10.1. PulseNet WGS Workflow



10.2. MiSeq Sequencing Workflow



LABORATORY STANDARD OPERATING PROCEDURE FOR WHOLE GENOME SEQUENCING ON MISEQ

Doc. No. PNL38

Ver. No. 03

Effective Date:

Page 17 of 31

11. RELATED DOCUMENTS:

Document Number	Title
PNL33	DNA Extraction and QC SOP
PNL34	Nextera XT Library Prep SOP
PNL34.W1	Nextera XT Library Prep Workbook
PNL35	Illumina DNA Prep (DNA Flex) SOP
PNL35.W1	Illumina DNA Prep Workbook, 96 CD/UD Indexes, separate sample sheets for MiSeq IEM, MiSeq LRM (vs 3 and earlier versions), MiniSeq LRM and iSeq LRM
PNL35.W2	Illumina DNA Prep Workbook, 24 CD Indexes, separate sample sheets for MiSeq IEM, MiSeq LRM (vs 3 and earlier versions) and iSeq LRM
PNL35.W3	Illumina DNA Prep Checklist
PNL36	QIAseq FX Library Prep SOP
PNL36.W1	QIAseq FX Library Prep Workbook
PNL37	KAPA HyperPlus Library Prep SOP
PNL37.W1	KAPA HyperPlus Library Prep Workbook
PNL42	NEB Next Ultra II FS Library Prep SOP
PNL42.W1	NEB Next Ultra II FS Library Prep Workbook
PNL43	Illumina DNA PCR-Free Prep SOP
PNL43.W1	Illumina DNA PCR-Free Workbook
PNL43.W2	Illumina DNA PCR-Free Checklist
PNQ07	Illumina Sequence Data QC SOP
PND20	BioNumerics RefID Database Workflow SOP
PND05	BioNumerics Organism-Specific Database Workflow SOP

12. REFERENCES:

- 12.1. Illumina, Inc. Illumina Experiment Manager Software Guide. (Doc.# 15031335 v08). December 2018. https://support.illumina.com/content/dam/illumina-support/documents/documentation/software_documentation/illumina-experiment-manager-software-guide-15031335-08.pdf
- 12.2. Illumina, Inc. MiSeq Reporter Software Guide. (Doc.# 15042295 v05). January 2017. https://support.illumina.com/content/dam/illumina-support/documents/documentation/software_documentation/miseqreporter/miseq-reporter-software-guide-15042295-05.pdf
- 12.3. Illumina, Inc. MiSeq System Guide. (Doc# 15027617v05). August 2019. https://support.illumina.com/content/dam/illumina-support/documents/documentation/system_documentation/miseq/miseq-system-guide-for-local-run-manager-15027617-05.pdf

LABORATORY STANDARD OPERATING PROCEDURE FOR WHOLE GENOME SEQUENCING ON MISEQ

Doc. No. PNL38

Ver. No. 03

Effective Date:

Page 18 of 31

- 12.4. Illumina, Inc. Illumina Experiment Manager: Narration Transcript.
https://support.illumina.com/content/dam/illumina-support/courses/illumina-experiment-manager/story_content/external_files/Experiment_Manager_Narration_Transcript.pdf
- 12.5. Illumina, Inc. Local Run Manager v2 Software Guide (Doc # 1000000002702 v08). November 2020. https://support.illumina.com/content/dam/illumina-support/documents/documentation/software_documentation/local-run-manager/local-run-manager-v2-software-guide-1000000002702-08.pdf
- 12.6. Illumina, Inc. Local Run Manager v3 Software Guide (Document # 1000000111492 v01). March 2022. https://support.illumina.com/content/dam/illumina-support/documents/documentation/software_documentation/local-run-manager/local-run-manager-v3-software-guide-1000000111492_01.pdf
- 12.7. Illumina, Inc. MiSeq System Custom Primers Guide (15041638 v01). March 2016. <https://emea.support.illumina.com/downloads/miseq-system-custom-primers-guide-15041638.html>
- 12.8. Illumina, Inc. Customer Notification. November 10, 2021. <https://support.illumina.com/bulletins/2019/09/impact-of-ammonium-based-cleaning-products-on-sequencing-run-per.html>
- 12.9. Illumina, Inc. How to edit a sample sheet and requeue an analysis in BaseSpace Sequence Hub (illumina.com). <https://support.illumina.com/bulletins/2018/03/how-to-edit-a-sample-sheet-and-requeue-a-miseq--hiseq--or-novase.html>

13. CONTACTS:

- 13.1. PulseNet NGS Lab:
pulsenetngslab@cdc.gov

14. AMENDMENTS:

- 14.1. **01/31/2019:**
 - Pulled procedures pertaining to setting up a run on the MiSeq out of PNL32 to create separate document (PNL38)
- 14.2. **02/02/2021**
 - Added additional definitions
 - Updated reagent and supply catalog numbers
 - Added information for the NEBNext Ultra II FS and QIAseq FX library prep
 - Added the “Review Run Metrics” section (9.8.)
 - Added the workflow charts
 - Added related documents
 - Updated references
 - Updated the PhiX instructions
 - Updated the MiSeq sample sheet instructions by adding LRM
- 14.3. **02/01/2023**
 - Added information for DNA PCR-Free Prep library prep kit
 - Updated Table 1 (thawing and post-thaw storage of the cartridge)

LABORATORY STANDARD OPERATING PROCEDURE FOR WHOLE GENOME SEQUENCING ON MISEQ			
Doc. No. PNL38	Ver. No. 03	Effective Date:	Page 19 of 31

- Additional figures added for clarity
- Added Table 7 summarizing the instrument maintenance
- Updated notes
- Updated references
- Former appendix PNL38-2 (MiSeq sample sheet set up) was separated into two appendixes by the software type: LRM (PNL38-2) and IEM (PNL38-3). Former appendix PNL38-3 (requeuing the analysis) is now appendix PNL38-4.
- Retired the MiSeq LRM job aid PNL38.JA1. This information is now included in the appendix PNL38-2.

15. APPROVAL SIGNATURES:

Approved By: _____ Date: _____
PulseNet QA/QC Personnel

Approved By: _____ Date: _____
PulseNet Outbreak Detection and Surveillance Unit Chief

Approved By: _____ Date: _____
PulseNet Next Generation Subtyping Methods Unit Chief

Approved By: _____ Date: _____
PulseNet Reference Outbreak Surveillance Team Lead

LABORATORY STANDARD OPERATING PROCEDURE FOR WHOLE GENOME SEQUENCING ON MISEQ			
Doc. No. PNL38	Ver. No. 03	Effective Date:	Page 20 of 31

16. APPENDICES:

Appendix PNL38-1

Addition of PhiX Control

1. Purpose

PhiX is a 500 bp balanced genome of a bacteriophage and is recommended by Illumina as an instrument control which may help indicate how well the sequencing instrument is performing. A PhiX run after instrument maintenance or repair is recommended.

2. Procedure

NOTE1: According to Illumina, 12.5 pM and 20 pM PhiX dilutions produce optimal cluster density for v2 and v3 reagents, respectively.

NOTE2: If using a previously prepared PhiX denatured dilution, proceed to Step 2.7.

- 2.1. Prepare 0.2 N NaOH by adding 400 µl of molecular-grade water to a 100 µl aliquot of 1 N NaOH and place on ice until ready to use.
- 2.2. Combine 1 µl of 10 nM PhiX library from the stock and 4 µl of RSB (or the diluent, e.g., EBT, Tris-HCl, your library prep kit uses to normalize libraries/library pool) to obtain 2 nM PhiX library in a new 1.5 ml microcentrifuge tube.
- 2.3. Combine 5 µl of 2 nM PhiX library and 5 µl of 0.2 N NaOH and mix well.
- 2.4. Incubate for 5 minutes at room temperature to denature the PhiX library into single strands.
- 2.5. **IMMEDIATELY** add 490 µl of pre-chilled HT1 to the tube containing 10 µl of denatured 1 nM PhiX library after the incubation to obtain 20 pM of denatured PhiX library.
- 2.6. Label, date and initial the microcentrifuge tube. Denatured 20 pM PhiX libraries may be stored up to three weeks at -15° to -25°C.
- 2.7. If using v2 sequencing reagents, dilute the denatured 20 pM PhiX library to 12.5 pM by adding 62.5 µl of 20 pM denatured PhiX library to 37.5 µl of pre-chilled HT1.
- 2.8. Determine percentage of PhiX (N%) to add to run. See below for general guidance.
 - If troubleshooting a library preparation: 1-10%
 - If troubleshooting a MiSeq instrument issue: 5% or above (depending on Illumina tech support guidance) is appropriate
 - If using PhiX routinely in every run: 1% is recommended
 - If loading libraries with low sample diversity, or spiking in amplicon: 5-10%
- 2.9. Combine N x 10µl of 12.5 pM (or 20 pM) PhiX with 1000 – (N x 10) µl of sample library of desired concentration.
 - For example: for 1% PhiX addition, add 10 µl of 12.5 pM PhiX to 990 µl of denatured 10-20 pM sample library.
- 2.10. Place on ice until ready to proceed with heat denaturation of sample library (see section 9.3.4.) and addition of PhiX.

3. PhiX Expected Results

- 3.1. The “Aligned (%)” in SAV/BaseSpace upon run completion should match the percentage of PhiX spiked into a run. If not, this may be an indication of an instrument or library issue.

LABORATORY STANDARD OPERATING PROCEDURE FOR WHOLE GENOME SEQUENCING ON MISEQ			
Doc. No. PNL38	Ver. No. 03	Effective Date:	Page 21 of 31

- 3.1.1. Libraries that are smaller than the size of the PhiX library (500 bp) will outcompete the PhiX, resulting in lower aligned percentage.

NOTE: *Libraries with low GC content, such as Campylobacter are prone to over-tagmentation particularly with the Nextera XT prep resulting in shorter fragments.*

- 3.1.2. Libraries that are larger than the PhiX will result in greater percentage of PhiX aligned than what was spiked in.

NOTE: *Libraries with high GC content, such as Mycobacterium tuberculosis tend to under-tagment particularly with the Nextera XT prep resulting in larger fragments.*

- 3.2. Contact Illumina Tech Support and cc pulsenetngslab@cdc.gov for troubleshooting assistance if necessary.

LABORATORY STANDARD OPERATING PROCEDURE FOR WHOLE GENOME SEQUENCING ON MISEQ

Doc. No. PNL38

Ver. No. 03

Effective Date:

Page 22 of 31

Appendix PNL38-2

MiSeq Sample Sheet Setup Using LRM

NOTE: *MiSeq sample sheet can be set up either by importing the SampleSheet tab from the library prep workbook (section 1) or by entering the information manually (section 2).*

1. Generate and Import a Sample Sheet from the PulseNet Workbook

- 1.1. Ensure that the Initial Dilution tab (Nextera XT) or the Library Prep tab (DNA Prep) is filled out properly. These tabs will auto populate the sample sheets with the necessary information. Check each field in the Sample Sheet tab to ensure that all of the data has been accurately populated based on the Initial Dilution (Nextera XT) or the Library Prep (Illumina DNA Prep) tab. The Description (Column J for XT or Column K for DNA Prep) and Sample_Project (Column I for XT or Column J for DNA Prep) fields may be blank.

NOTE: *The index sequences are auto populated from the Indices tab. Do not modify or delete this tab.*

- 1.2. Convert the appropriate Sample Sheet tab to a .csv (either by saving as or exporting as a .csv comma delimited file)

NOTE: *The appropriate sample sheet tab is determined by the MiSeq Control Software installed on the MiSeq. See Table 8 below for which sample sheet tab to use based on the version of MCS in use.*

MCS Version	Run Set Up Software	Sample Sheet Tab to Use
v4*	Local Run Manager (v 3.0 and above)	LRM3 SampleSheet
v3	Local Run Manager (v 2.0)	LRM SampleSheet
2.6.2 and below	IEM	IEM_SampleSheet – see Appendix PNL 38-3 for instructions
*Instruments upgraded to Windows 10 will have MCS v4 installed and the LRM3 SampleSheet tab should be used.		

Table 8. Version of MCS, run set up software and corresponding sample sheet tabs to be used for run set up on the MiSeq.

1.2.1. To save as .csv:

- 1.2.1.1. Open the “File” menu option and select “Save As”.
- 1.2.1.2. Change “Save as type:” to “CSV (Comma delimited)”.
- 1.2.1.3. Navigate to the desired file folder location and save the file using the Plate ID as the filename (i.e., LabID-MXXXX-YYMMDD).
- 1.2.1.4. Click “OK” to save only the active sheet, and “Yes” in the following window to keep using the CSV format.

OR

1.2.2. To export as .csv:

- 1.2.2.1. Open the “File” menu option and select “Export”.
- 1.2.2.2. Select “Change File Type”.

LABORATORY STANDARD OPERATING PROCEDURE FOR WHOLE GENOME SEQUENCING ON MISEQ

Doc. No. PNL38

Ver. No. 03

Effective Date:

Page 23 of 31

- 1.2.2.3. Choose “CSV (Comma delimited)”.
- 1.2.2.4. Save and select “OK” to save only the active sheet.
- 1.3. Open the CSV file in WordPad or Notepad and delete any commas listed after the final sample. See Fig. 6 below.

NOTE: *This must be done in order to avoid errors when importing the sample sheet.*




Figure 6. Sample sheet tab in Notepad, depicting extra commas to be deleted.

- 1.4. Transfer the .csv file to the instrument using a flash drive or similar.
 - 1.4.1. On the MiSeq, open Chromium.
 - 1.4.2. On the LRM home page, select “Create Run”.
 - 1.4.3. Select the “GenerateFASTQ” analysis module (Fig. 7).




Figure 7. GenerateFASTQ module option in Local Run Manager.

- 1.4.4. Select “Import Sample Sheet” and navigate to the .csv file generated from the workbook.
- 1.4.5. Ensure that ALL run parameters, sample and index information imported correctly (Fig. 8). This includes Run ID, Library Prep Kit, Index Kit, **Adapter Trimming**, Read Length, **Custom Primers (if sequencing Illumina DNA PCR-Free Prep libraries)**, etc.

NOTE: *If the Library Prep Kit defaults to “Custom” (as in the example in Fig. 9 below); do not proceed. Ensure that the commas in the exported sample sheet were deleted (Step 1.3)*

LABORATORY STANDARD OPERATING PROCEDURE FOR WHOLE GENOME SEQUENCING ON MISEQ

Doc. No. PNL38

Ver. No. 03

Effective Date:

Page 24 of 31

	SAMPLE ID *	SAMPLE DESCRIPTION	PLATE *	INDEX WELL *	INDEX 1 (I7)	INDEX 2 (I5)	SAMPLE PROJECT	
1	03-98-MS-M3235-22-038		D	F11	UDP0374	UDP0374		x
2	C6472-MS-M3235-22-038		D	G11	UDP0375	UDP0375		x
3	H89394-MS-M3235-22-038		D	H11	UDP0376	UDP0376		x
4	D5480-MS-M3235-22-038		D	A12	UDP0377	UDP0377		x

Figure 8. Screen shot of LRM with run parameters correctly populated for Illumina DNA Prep 500c run.

	SAMPLE ID *	SAMPLE DESCRIPTION	PLATE *	INDEX WELL *	INDEX 1 (I7)	INDEX 2 (I5)	SAMPLE PROJECT	
1	03-98-MS-M3235-22-038		D	F11	UDP0374	UDP0374		x
2	C6472-MS-M3235-22-038		D	G11	UDP0375	UDP0375		x
3	H89394-MS-M3235-22-038		D	H11	UDP0376	UDP0376		x
4	D5480-MS-M3235-22-038		D	A12	UDP0377	UDP0377		x

Figure 9. Screen shot of Local Run Manager with Library Prep kit type defaulting to “Custom”.

- 1.4.6. Upon confirming all of the run parameters and sample information is correct, select “Save Run”. The run will now be displayed on the LRM home page.
- 1.4.7. Close Chromium. The run will now be saved in LRM and will show up on the MiSeq Control Software homepage for sequencing.

LABORATORY STANDARD OPERATING PROCEDURE FOR WHOLE GENOME SEQUENCING ON MISEQ			
Doc. No. PNL38	Ver. No. 03	Effective Date:	Page 25 of 31

2. Manually Enter Run Information in LRM:

2.1. On the MiSeq, open Chromium.

2.2. On the LRM home page, select “Create Run” to set up run parameters.

2.3. Select the “GenerateFASTQ” analysis module.

2.4. Manually enter run information:

2.4.1. Run Name: enter the Plate ID (i.e., LabID-MXXXX-YYMMDD).

2.4.2. Run Description

2.4.3. Run Settings:

2.4.3.1. Library Prep Kit: select appropriate kit in drop-down.

2.4.3.2. Read Type: select “Paired End”.

2.4.3.3. Index Reads: select “2”.

2.4.3.4. Read lengths: input desired number of cycles, depending on kit. For 500 cycles, type “251” for READ 1 and READ 2. For 300 cycles, type “151” for READ 1 and READ 2.

2.4.3.5. **If sequencing Illumina DNA PCR-Free libraries:** Custom primers: “CustomRead1PrimerMix (C1)”.

2.4.4. **Adapter Trimming:** confirm box is checked.

2.4.5. Input sample information, including sample ID and indexes used.

2.5. Select “Save Run”. The run will now be displayed on the LRM home page.

2.6. Close Chromium.

LABORATORY STANDARD OPERATING PROCEDURE FOR WHOLE GENOME SEQUENCING ON MISEQ			
Doc. No. PNL38	Ver. No. 03	Effective Date:	Page 26 of 31

Appendix PNL38-3

MiSeq Sample Sheet Setup using IEM

NOTE: *MiSeq sample sheet can be set up either by importing the SampleSheet tab from the library prep workbook (section 1) or by entering the information manually into IEM (section 2).*

1. Generate and Import a Sample Sheet from the PulseNet Workbook

- 1.1. Ensure that the Initial Dilution tab (Nextera XT) or the Library Prep tab (DNA Prep) is filled out properly. These tabs will auto populate the sample sheets with the necessary information. Check each field in the Sample Sheet tab to ensure that all of the data has been accurately populated based on the Initial Dilution (Nextera XT) or the Library Prep (Illumina DNA Prep) tab. The Description (Column J for XT or Column K for DNA Prep) and Sample_Project (Column I for XT or Column J for DNA Prep) fields may be blank.

NOTE: *The index sequences are auto populated from the Indices tab. Do not modify or delete this tab.*

- 1.2. Convert the appropriate Sample Sheet tab to a .csv (either by saving as or exporting)

1.2.1. To save as .csv:

- 1.2.1.1. Open the “File” menu option and select “Save As”.
- 1.2.1.2. Change “Save as type:” to “CSV (Comma delimited)”.
- 1.2.1.3. Navigate to the desired file folder location and save the file using the Plate ID as the filename (i.e., LabID-MXXXX-YYMMDD).
- 1.2.1.4. Click “OK” to save only the active sheet, and “Yes” in the following window to keep using the CSV format.

OR

1.2.2. To export as .csv:

- 1.2.2.1. Open the “File” menu option and select “Export”.
- 1.2.2.2. Select “Change File Type”.
- 1.2.2.3. Choose “CSV (Comma delimited)”.
- 1.2.2.4. Save and select “OK” to save only the active sheet.

- 1.3. Open the CSV file in WordPad or NotePad and delete any commas listed after the final sample.

NOTE: *This must be done in order to avoid errors when importing the sample sheet.*

- 1.4. Transfer the .csv file to the instrument using a flash drive or similar.

- 1.4.1. Open the Illumina Experiment Manager (IEM) software, select “Edit Sample Sheet” and select the CSV file generated from the workbook.

- 1.4.2. In the new window, copy and paste the CSV file among the list of other sample sheets and select “Open”.

- 1.4.3. In the “Run Settings” window, verify that all information is correct.

- 1.4.3.1. Enter the Plate ID as the “Experiment Name” and check the box for “**Use adapter trimming**”.

- 1.4.3.2. **If sequencing Illumina DNA PCR-Free Prep libraries**, check the box for “CustomRead1PrimerMix” (C1) (right side of screen).

- 1.4.3.3. Click “Next”.

LABORATORY STANDARD OPERATING PROCEDURE FOR WHOLE GENOME SEQUENCING ON MISEQ

Doc. No. PNL38

Ver. No. 03

Effective Date:

Page 27 of 31

- 1.4.4. Verify the Sample Sheet is valid by checking for the green “Valid” status in the “Sample Sheet Status” field on the left-hand side of the screen.
- 1.4.5. Click “Finish”.
- 1.4.6. Save the Sample Sheet using the same Plate ID (i.e., LabID-MXXXX-YYMMDD).
- 1.4.7. Click “Yes” to overwrite the previous file and “No” when prompted to view the sheet in Excel.

2. Manually Enter Run Information in IEM:

- 2.1. Create MiSeq Sample Plate – The Sample Plate stores information in a 96-well plate format, designates index assignment and allows creation of a Sample Sheet in the proper format.
 - 2.1.1. Open Illumina Experiment Manager (IEM icon on the MiSeq screen).
 - 2.1.2. Choose “Create Sample Plate”.
 - 2.1.3. Select the appropriate Sample Prep Kit and click “Next” (e.g., NexteraXT).
 - 2.1.4. Enter the Plate ID (i.e., LabID-MXXXX-YYMMDD).
 - 2.1.5. Select “2 (Dual)” for the number of index reads and select “Next”.
 - 2.1.6. Under the “Plate” tab, start in position A01 and enter the first Sample ID (PulseNet BioNumerics Key). Repeat for all samples.
 - 2.1.7. Select “Sample Name” and enter the Sample Name which should be the Sample ID plus the Plate ID (<Sample ID> - <Plate ID>). Repeat for all samples.
 - 2.1.8. Assign indexes for all samples by selecting “Index Well” or “Index1 (I7)” and “Index2 (I5)” from the drop-down menu at the top left of the screen. Select the appropriate index from the drop-down menu in each cell for all samples.
 - 2.1.9. Verify the layout by clicking on the “Table” tab, confirming that information in the wells correspond to their positions on the library prep worksheet.
 - 2.1.10. Click “Finish”.
 - 2.1.11. Save the Plate file as the Plate ID (i.e., LabID-MXXXX-YYMMDD).
- 2.2. Create the Sample Sheet – The Sample Sheet sets the parameters for run set-up, performance and analysis.
 - 2.2.1. Open the Illumina Experiment Manager (IEM icon on the MiSeq screen).
 - 2.2.2. Choose “Create Sample Sheet”.
 - 2.2.3. Choose “MiSeq” and select “Next”.
 - 2.2.4. Select the appropriate category and application and then press the “Next” button.
 - 2.2.4.1. Category: Other.
 - 2.2.4.2. Application: FastQ Only
 - 2.2.5. On the “Run Settings” page, enter the following:
 - 2.2.5.1. Reagent Cartridge Barcode: provide Plate ID (i.e., LabID-MXXXX-YYMMDD).
 - 2.2.5.2. Library Prep Workflow: Select the library prep method used (e.g., Nextera XT).
 - 2.2.5.3. Index Adapters: Select the indices used for library prep.
 - 2.2.5.4. Index Reads: 2 (Dual).
 - 2.2.5.5. Experiment Name: Enter the Plate ID.
 - 2.2.5.6. Investigator Name: Initials of technician(s) performing the run.
 - 2.2.5.7. Description: May use the Plate ID.
 - 2.2.5.8. Cycle Reads:
 - 2.2.5.8.1. For a v2 300 cycle kit, set at 151 x 151.

LABORATORY STANDARD OPERATING PROCEDURE FOR WHOLE GENOME SEQUENCING ON MISEQ			
Doc. No. PNL38	Ver. No. 03	Effective Date:	Page 28 of 31

- 2.2.5.8.2. For a v2 500 cycle kit, set at 251 x 251.
- 2.2.5.8.3. For a v3 600 cycle kit, set at 251 x 251 (recommended) or 301 x 301.
- 2.2.5.9. Leave all other settings as default.
- 2.2.5.10. Verify that the “**use adapter trimming**” is checked
- 2.2.5.11. **If sequencing Illumina DNA PCR-Free Prep libraries**, ensure “CustomRead1PrimerMix” (C1) box is checked (right-hand side of screen).
- 2.2.5.12. Select “Next”.
- 2.2.6. On the “Sample Selection” page, click on “Select Plate”.
NOTE: *If this option is not displayed, uncheck the “Maximize” box at the top right-hand side of the screen.*
- 2.2.7. Find the appropriate plate from the directory in the new pop-up window.
- 2.2.8. Highlight the correct plate and select “Open”.
- 2.2.9. Click the “Select All” button at the bottom left-hand side of screen.
- 2.2.10. Click “Add Selected Samples”.
- 2.2.11. Verify the successful addition of all desired samples to the Sample Sheet in the “Samples To Include in the Sample Sheet Field” window on the right-hand side of the screen.
- 2.2.12. Verify the Sample Sheet is ‘valid’ by checking for the green “Valid” status in the “Sample Sheet Status” field on the right-hand side of the screen.
- 2.2.13. Click “Finish”.
- 2.2.14. Save the Sample Sheet using the Plate ID (i.e., LabID-MXXXX-YYMMDD). Select “No” when prompted to view the sheet in Excel.

LABORATORY STANDARD OPERATING PROCEDURE FOR WHOLE GENOME SEQUENCING ON THE MISEQ			
Doc. No. PNL38	Ver. No. 02	Effective Date:	Page 29 of 31

Appendix PNL38-4

Requeuing a Run for Analysis

1. Purpose

Requeuing a run may be necessary if an indexing or sample naming error has been identified (i.e., the index assignments in the sample sheet are not accurate). Prior to requeuing analysis, prepare a corrected sample sheet. This will be used by the instrument for the new analysis, to re-parse the data and to rename the associated fastq files. A run can be requeued for analysis either using MiSeq Reporter (section 2.1.) or LRM (section 2.2.) or BaseSpace (section 2.3.).

2. Procedure

2.1. Requeue a Run Using MiSeq Reporter (MSR)

- 2.1.1. In order to apply the corrected or new analysis parameters to a run, a new sample sheet must first be created. To select a run to be re-queued with a new sample sheet, it is important that the input files are present in the analysis run folder. SampleSheet.csv, RTAComplete.txt, and RunInfo.xml must all be present.
- 2.1.2. In the analysis run folder (D:\ Data\Illumina\MiSeq Analysis/"Recent run") rename the old sample sheet (for example, "SampleSheet.csv.old") and paste the newly revised sample sheet in the analysis folder under the name "SampleSheet.csv".
 - 2.1.2.1. If correcting index assignment, the new sample sheet should have the corrected index pair assignments and corresponding index sequences.
- 2.1.3. Connect to MiSeq Reporter on the local computer by opening Internet Explorer. If MiSeq Reporter does not open with the launch of Internet Explorer, type the following in the browser bar: <http://localhost:8042/>.
- 2.1.4. Check the checkbox next to the run and Click "Requeue". The analysis will take roughly an hour depending on the extent of changes to the Sample Sheet.

NOTE1: *Alternatively, the "QueueForAnalysis.txt" file can be deleted in the MiSeq Analysis/"recent run folder" to begin the reanalysis of the dataset.*

NOTE2: *MiSeq Reporter requires the following files which are generated during the sequencing run to perform secondary analysis or to requeue analysis: all *.bcl, *.filter, and *.locs.*
- 2.1.5. The reanalyzed data (i.e., fastq files) will be stored in the "Basecalls" folder of the run analysis and output folders.
- 2.1.6. Additional detailed instructions on how to requeue analyses from the MSR interface may be found in the MiSeq Reporter Software Guide, page 11:
https://support.illumina.com/content/dam/illumina-support/documents/documentation/software_documentation/miseqreporter/miseq-reporter-software-guide-15042295-05.pdf

2.2. Requeue a Run Using LRM and MCS v3.0 or More Recent

- 2.2.1. Open Chromium to access LRM.
- 2.2.2. Find the run to be requeued.
- 2.2.3. Click "Actions".

LABORATORY STANDARD OPERATING PROCEDURE FOR WHOLE GENOME SEQUENCING ON THE MISEQ

Doc. No. PNL38

Ver. No. 02

Effective Date:

Page 30 of 31

2.2.4. Click “Requeue”.

2.2.5. Select “Edit Setup” when prompted by the software.

2.2.6. Import the new sample sheet and proceed.

2.2.7. The reanalyzed data (i.e., fastq files) will be stored in the “Alignment_2” folder of the run analysis and output folders.

NOTE: *If a run was accidentally completed with the Library Prep Kit listed as “Custom”, indexes and adapter trimming cannot be edited using requeue for analysis as described above. Instead, create a new run in LRM (recommended to use the same Run ID as the run that needs requeuing, with a “v2” at the end) and import the corrected sample sheet (with the new “v2” Run ID). After setting up this run in LRM, right click and select “Import” on the LRM home screen (Fig.10). Copy the file directory for the run that needs to be re-analyzed into the dialog box (Fig. 11); ensure that the output folder is correct and select “Import Run”. The resulting run will have properly demultiplexed and adapter trimmed fastq files.*




Figure 10. The drop-down menu to select “Import” for a manual re-analysis of run data.




Figure 11. Selecting the run to be re-analyzed and the output folder in the “Import run” pop-up window.

LABORATORY STANDARD OPERATING PROCEDURE FOR WHOLE GENOME SEQUENCING ON THE MISEQ			
Doc. No. PNL38	Ver. No. 02	Effective Date:	Page 31 of 31

2.3. Requeue a Run Using BaseSpace Sequence Hub

NOTE1: *In order to requeue the analysis within BaseSpace, you must be the owner of the run. A shared run does not allow permissions to requeue analysis. Additionally, **a requeued run can only be submitted up to five times in BaseSpace Sequence Hub**. If assistance with additional requeues is needed, email Illumina Technical Support, and be prepared to share your run ID and sample sheet.*

NOTE2: *MSR and BaseSpace will requeue a new sample sheet independent of each other. If a new sample sheet is requeued with MSR, it will not update the information on BaseSpace.*

2.3.1. Follow the general instructions for requeuing in BaseSpace on the Illumina website:

<https://support.illumina.com/bulletins/2018/03/how-to-edit-a-sample-sheet-and-requeue-a-miseq--hiseq--or-novase.html>

2.3.1.1. In order to fix the **sample sheet**, refer to the following instructions at the Illumina webpage: <https://help.basespace.illumina.com/runs/fix-sample-sheet>

2.3.1.2. In order to fix the **index assignment**, refer to the following instructions at the Illumina webpage: <https://help.basespace.illumina.com/sequence/fix-indexes>