




Center for Community and Preventive Health
Infectious Disease Epidemiology Section
Influenza Virologic Surveillance Handbook
2013 – 2014 Season

Version 1.0



September 2013

Mutation

Antigenic drift




Antigenic shift



I. INTRODUCTION

Virologic surveillance is the foundation on which national and international influenza surveillance systems are built. The goal of virologic surveillance is to identify and track drift variants of currently circulating influenza virus types and subtypes and to detect the emergence of novel influenza A subtypes in human populations. This information allows for monitoring of the match between vaccine strains and currently circulating viruses and selection of optimal vaccine components each year.

II. DATA FLOW OF THE U.S. VIROLOGIC SURVEILLANCE SYSTEM



III. LOUISIANA VIROLOGIC SURVEILLANCE

Beginning in the 2013-2014 influenza season, the goal is for the Louisiana State Public Health Laboratory to increase samples to meet requirements of the Association of Public Health Laboratories *Influenza Virologic Surveillance Right Size Roadmap*. The increase in sample submission will require regular participation from a core group of surveillance sites statewide. All materials required for sample collection and submission will be provided free of charge and transportation will be coordinated through Statewide transport.

Participation in active surveillance will require:

- Collecting a nasal or nasopharyngeal (NP) swab on **all patients** who present with clinical symptoms resembling influenza-like illness on any one day of the week (or more if the site is willing).
- Packing specimens in an ice chest with proper ice blocks (all provided) for Statewide transport pick-up.

A portion of flu positives from active surveillance will be forwarded to CDC for further antigenic characterization and antiviral resistance testing. All **flu negative NP swabs** submitted will be tested for

other respiratory viruses including Respiratory Syncytial Virus, Parainfluenza, Human Metapneumovirus, and Adenovirus.

A. Specimen Collection


Nasopharyngeal (NP) swabs, when collected by properly trained personnel are the optimal choice of clinical material. A nasal swab is an acceptable substitute and is the optimal choice for many participating sites.

Nasopharyngeal (NP) swabs

- 1) Optimal timing. Specimens ideally should be collected within 72 hours of influenza-like illness symptom onset (e.g. respiratory symptoms and/or fever) but are acceptable up to 5 days from symptom onset. Specimens should ideally be collected prior to the initiation of antiviral medications but are acceptable after the antiviral therapy has begun.
- 2) Materials. Nasopharyngeal swab (flexible shaft) with rayon tip, viral transport medium. All materials will be provided for participants.
- 3) Swab types. Swab specimens should be collected using only swabs with a synthetic tip, such as nylon or Dacron®, and an aluminum or plastic shaft. Calcium alginate swabs are unacceptable and cotton swabs with wooden shafts are not recommended.
- 4) Collecting the NP swab. Insert swab through the nares parallel to the palate (not upwards, Figure 1) until resistance is encountered or the distance is equivalent to that from the ear to the nostril of the patient. An instructional video from The Joint Commission on NP swab collection can be viewed at:

<http://www.youtube.com/watch?v=hXohAo1d6tk&feature=youtu.be>

Figure 1: Nasopharyngeal Swab Collection



B. Specimen Packaging

Refrigerate specimens after collection. Ship specimens refrigerated to be received within 72 hours from collection.




Use a collection kit containing swabs and viral transport media to collect flu specimen.



The swab on the left is for nasal specimens. The swab on the right is for nasopharyngeal.



After collection, place the swab in the tube of viral transport media.



Snap the swab shaft and cap the tube. Double check that your media is in date.



Put the date and time of collection on the tube along with the patient's name.



For each specimen collected, completely fill out a Lab Test Request Form 96.



Place each specimen in a separate compartment of the bubble wrap pouch.



Roll up the bubble wrap sleeve with the specimens. Tape the roll closed.



Each ice chest is furnished with two biohazard transport bags and absorbent strips.



Insert specimen roll into the long pouch. Pull off the liner and press to close.



Fold the lab forms in half and insert them into the outside pouch of the transport bag.



Freeze the ice bricks prior to use. Cover the bottom of the ice chest with one brick.



Place the transport bag with samples on top of the first ice brick.



Layer the second ice brick on top of the transport bag with samples.



Tape the ice chest closed. Ship ice chests overnight via StateWide Transport.


C. Lab Test Request Form

Lab Form 96 – Virology needs to accompany all specimens submitted for testing. When you enroll as a participating site, you will be provided with a form with your site information pre-filled.

Lab Form 96 Revision 12/2009		Lab Test Request Form Virology	
(*) INDICATES REQUIRED INFORMATION. INCOMPLETE INFORMATION MAY CAUSE SPECIMEN REJECTION.			
Patient Information		Submitter Information	
*First Name _____	*Last Name _____	*Facility Name _____	
*Date of Birth _____	*Gender _____	*Facility Address _____	
Patient's Home Address _____		*City, State _____	*Zip Code _____
City, State _____	Parish _____	Name of Contact Person _____	
Medicaid Number _____	Patient ID Number _____	Phone (____) _____	Fax (____) _____
Specimen Information			
Submitted Specimen Is From: ☒ Human ☐ Animal _____			
*Collection Date: _____ Time: _____		If Frozen, Indicate Date: _____ Time: _____	
*Specimen Type: ☒ Swab ☐ Aspirate/Wash ☐ Tissue ☐ Viral Culture			
*Specimen Source: ☐ Nasal ☐ Nasopharynx ☐ CSF ☐ Blood			
☐ Oropharynx ☐ Bronchil ☐ Acute Serum ☐ Stool			
☐ Other ☐ Trachea ☐ Convalescent Serum ☐ Vomitus			
Additional Information: ☐ Mother ☐ Follow-Up ☐ Prenatal ☐ Date of Symptom Onset: _____			
☐ Child ☐ Prenatal			
Test Request Information			
All tests listed may not currently be available. For questions regarding test availability, contact the Virology Dept. at 504-219-4676.			
*At Least One Test Must Be Requested			
☒ Respiratory Virus Panel* RSV, Influenza, Parainfluenza Metapneumovirus, Rhinovirus and Adenovirus		☐ Arbovirus Panel, IFA SLE IgM, SLE IgG, EEE, IgM, EEE IgG, CE IgM, CE IgG, WEE IgM and WEE IgG	
☐ Norovirus Real Time RT-PCR Norovirus GI and GII		☐ Arbovirus Panel, MIA West Nile and SLE	
☒ Influenza Real Time RT-PCR Detection and Characterization		☐ Maternal Serum Panel HCG, uE3, AFP and Inhibin A	
Epi Risk Factor <u>Influenza Piglet</u> <u>Size Viralologic Survey Marker,</u> <u>Sentinel Provider</u>		☐ Hepatitis A Hepatitis A Total Antibody (Anti-HAV) Hepatitis A IgM Antibody (Anti-HAV IgM)	
		☐ Hepatitis B Panel Hepatitis B Surface Antigen (HBsAg) Hepatitis B Core Total Antibody (Anti-HBc)	
		☐ Hepatitis B Immunization (Anti-HBs)	
		☐ Hepatitis B Core IgM Antibody (Anti-HBc IgM)	
		☐ Hepatitis C Total Antibody (Anti-HCV)	
		☐ Rubella IgG ☐ Mumps IgG	
		☐ Toxo IgG ☐ Vaccinia	
		☐ CMV IgG ☐ Varicella	
		☐ Measles IgM ☐ Other _____	
		☐ Measles IgG	
		☐ Other Testing ☐ Herpes IgM ☐ Herpes IgG ☐ Lyme IgM ☐ Lyme IgG ☐ Lyme Total Ab	
To Be Forwarded to CDC: Contact 504-219-4646 for prior approval.		☐ Yes ☐ No	
		Diagnosis Suspected: _____ Please include One CDC History form per specimen.	
Send To: DHH-OPH Central Lab, 3101 West Napoleon Ave., Metairie, LA 70001			
TO BE COMPLETED BY STATE LABORATORY			
LABORATORY NUMBER: _____		DATE/TIME RECEIVED: _____	
		TEMPERATURE CONDITION: _____	

D. Shipping Specimens


Website: **statewidetransportation.com**



Enter User ID: **centrallab**

Enter Password: **central**

Press **LOGIN**



Press ORDER and the ordering screen will open

Order - Windows Internet Explorer

http://107.7.3.210/cxt/order.asp

Order

Statewide Transport Inc

Logout

Order

Track

Account

Address Book

Required fields display asterisk(*)

Origin Find Origin Address

Country*: United States

Name*:

Address*:

Suite:

City*: State/Province*:

Postal Code*: -

Phone:

Remarks: (You may enter up to 900 characters) 900 characters left

Destination Find Destination Address

Country*: United States

Name*:

Address*:

Suite:

City*: State/Province*:

Postal Code*: -

Attention To:

Phone:

Done

Internet | Protected Mode: Off

75%

Press **FIND ORIGIN ADDRESS** and a screen will appear asking for the first few characters of the address you wish to find. Type the first few characters of your address and press **FIND ADDRESS**.

The address will be displayed. Press **SELECT**. This will auto fill the address field with your address.

Press **FIND DESTINATION ADDRESS** and a screen will appear asking for the first few characters of the address you wish to find. Type **LOU** and press **FIND ADDRESS** for the Central Lab address.

Press **PARCELS**.

http://107.7.3.210/cxt/switchaccount.asp?ShowParcelForm=list - Windows Internet Explo...

http://107.7.3.210/cxt/switchaccount.asp?ShowParcelForm=list

After you have finished adding parcels, you can click this window's close button (the x in the top right corner) and your selected parcels will automatically be saved.

Parcel List				
ID	Type	Weight	Barcode	Comment
	Please choose a type ▼			

[Add](#)

Done Internet | Protected Mode: Off 75%

Press ADD. If you have more than one parcel going to the same address, press ADD for every parcel. You can then close this screen.

You do not need to worry about weight or type or billing. You also do not need to worry about a pickup date and time if you are printing the labels when you are ready to ship.

Press CONTINUE

Do not choose the RETURN option. RETURN is for Central Lab to send back empty ice chests.

Print the waybills. You will get one for each parcel (ice chest)

Internet Shipping: View/Print Label - W...
http://107.7.3.210/cxt/shippingLabel.asp



DHH CENTRAL LAB
DHH-OPH CENTRAL LABORATORY
3101 W NAPOLEON AVE
METAIRIE LA 70001
ORDER: 556360
2013-03-06 14:41

Parcel Number: 2691177
Parcel Type:
Weight: 1
L x W x H: 0 x 0 x 0
Reference: X-0160-B2691177

SHIP TO:
LAFOURCHE HEALTH UNIT
2535 VETERANS BLVD
THIBODAUX LA 70301


X-0160-B2691177

PARCEL TRACKING NUMBER: 2691177

Internet | Protected Mode: Off 75%

and one return (for the driver to sign)

Internet Shipping: View/Print Label - W...
http://107.7.3.210/cxt/shippingLabel.asp



DHH CENTRAL LAB
DHH-OPH CENTRAL LABORATORY
3101 W NAPOLEON AVE
METAIRIE LA 70001

Return
2013-03-06 14:41
Pieces: 1
Weight: 1

SHIP TO:
LAFOURCHE HEALTH UNIT
2535 VETERANS BLVD
THIBODAUX LA 70301


270301

TRACKING NUMBER: 556360


X-0160-0556360

References: _____
Signature: _____

Internet | Protected Mode: Off 75%

The Statewide IT contact is Lambert. His cell numbers are 504-416-5158 and 985-662-4325.

If you are not already in the State Public Health Laboratory StarLIMS database, you will be asked to fill out a secure fax form.

Page 11 of 11