

From Pipes to Patterns: Tracking Antimicrobial Resistance in Hawaiian Wastewater



#IDLabCon

Doris Di, Ph.D.

Lead Wastewater Microbiologist

State Laboratories Division (SLD)

Hawaii State Department of Health (HDOH)



DISCLOSURE

No conflict of interest to declare

Content Page



1.Introduction

2.Pilot Study: Experiment design and workflow

3.Results & Discussion

4.Summary & Conclusion


Introduction



#IDLabCon

- What is wastewater surveillance?
- Wastewater Surveillance Program in Hawaii
- Why tracking Antimicrobial Resistance?

What is Wastewater Surveillance?



History of Wastewater Surveillance




1950s



Salmonella

Enteric bacterial
pathogen
outbreaks


2020




SARS-CoV-2

Pandemic
respiratory
virus

Today



SARS-CoV-2,
Influenza, RSV



Hepatitis A




Candida auris



Norovirus,
Rotavirus



Fentanyl, Heroin, Cocaine,
methamphetamine




PFAS

Bacteria, Virus,
Yeast, Illicit Drugs,
Chemicals....


#IDLabCon




Wastewater Surveillance Program in Hawaii



Wastewater Surveillance Program in Hawaii



City & County



9 Wastewater Treatment Plants
(serving **68.7%** of the total
population in the State of Hawai'i)

Wastewater Surveillance Data Report




Published on DOH Website Biweekly



Link:

<https://health.hawaii.gov/coronavirusdisease2019/what-you-should-know/covid-19-data-reports/>




#IDLabCon



Why tracking Antimicrobial Resistance?



Antibiotic resistance threats in the US



~50 000 deaths attributed to **Carbapenem-resistant** *Acinetobacter* and *Enterobacterales* globally

Antibiotic resistance threats list

Urgent threats

- **Carbapenem-resistant *Acinetobacter***
- *Clostridioides difficile*
- **Carbapenem-resistant *Enterobacterales***
- Drug-resistant *Neisseria gonorrhoeae*

Serious threats

- Drug-resistant *Campylobacter*
- Drug-resistant *Candida*
- ESBL-producing *Enterobacterales*
- Vancomycin-resistant *Enterococci* (VRE)
- Multidrug-resistant *Pseudomonas aeruginosa*
- Drug-resistant non-typhoidal *Salmonella*
- Drug-resistant *Salmonella* serotype Typhi
- Drug-resistant *Shigella*
- Methicillin-resistant *Staphylococcus aureus* (MRSA)
- Drug-resistant *Streptococcus pneumoniae*


Concerning threats

- Erythromycin-resistant group A *Streptococcus*
- Clindamycin-resistant group B *Streptococcus*

Watch-list

- Drug-resistant *Mycoplasma genitalium*
- Drug-resistant *Bordetella pertussis*

Why tracking Antimicrobial Resistance?



In Hawaii, clinical cases increased in recent years, especially KPC and NDM.

Source: Disease Outbreak Control Division (DOCD), HDOH

#IDLabCon



Pilot Study: Experiment Design and Workflow



#IDLabCon

Study Sites and Sampling Time



- Pilot Location: **Hawaii Island and Kauai**



- WWTP 1: population served ~32,000
- WWTP 2: population served ~10,000

- Study Period:


- Jan 1st to March 31st, 2024 (13 weeks)
- 1 sample/week (duplicate)

- Objective:

- To understand the spatial-temporal distribution and genetic diversity of AMR genes in the two Hawaiian communities.



Workflow and Methods





**35 mL (×2)
Wastewater
Samples**



**Ceres Nanotrap®
Microbiome B Particles
+
Enhancement Reagent 3**



**QIAGEN AllPrep®
PowerViral®
DNA/RNA Kit**




**QIAGEN CLC Genomics
Workbench**



**QIAseq® FX DNA Library Kit
QIAseq® xHYB AMR Panel (2,786 ARGs)
&
Illumina MiSeq™ System
MiSeq™ Reagent Kit v2
Illumina NextSeq 2000 P1 (300 cycles)**



**QIAcuity® Digital PCR 26k 24-well
&
dPCR Microbial DNA Detection Assays
(NDM, IMP, VIM, KPC, and OXA-48)**




Results & Discussion





#IDLabCon

- Temporal and Spatial Trends
- AMR Metagenomic Diversity
- Bacterial Diversity

Results: Temporal and Spatial Trends



Results: Temporal and Spatial Trends




Over a 13-Week period, dPCR showed that

- **Both WWTPs** showed **high concentration** of **KPC** and **OXA-48** genes.
 - WWTP 1: KPC ($7.0 \times 10^3 - 3.8 \times 10^4$ GC/mL)
 - WWTP 1: OXA-48 ($7.8 \times 10^2 - 4.0 \times 10^4$ GC/mL)
 - WWTP 2: KPC ($1.9 \times 10^4 - 5.7 \times 10^4$ GC/mL)
 - WWTP 2: OXA-48 ($1.6 \times 10^1 - 1.4 \times 10^2$ GC/mL)
- **IMP** genes were **frequently being detected** in **both WWTPs** but at a **low concentration** (1.2 – 68.9 GC/mL and 1.0 – 1.9 GC/mL).
- **NDM** genes were **consistently being detected** in **WWTP 1** (2.3 – 15.8 GC/mL) but were only detected two times in WWTP 2 (1.4 & 1.7 GC/mL).
- **VIM** genes were **detected more frequently** in **WWTP 2** (1.2 – 4.7 GC/mL) than in WWTP 1 (1.4 GC/mL) at low concentration.




Results: AMR Metagenomic Diversity



WWTP 1

WWTP 2




#IDLabCon



Results: AMR Metagenomic Diversity



WWTP 1



KPC ✓

OXA-48 ✓

IMP ✓

NDM ✓


VIM ✗

WWTP 1							
Carbapenem Resistance Genes of Interest [Abundance (%) Rank]							
Week 1		Week 5		Week 9		Week 13	
OXA beta-lactamase	7.78	OXA beta-lactamase	4.92	OXA beta-lactamase	3.07	OXA beta-lactamase	4.66
KPC beta-lactamase	1.91	OXA-12	1.98	KPC beta-lactamase	2.03	KPC beta-lactamase	1.78
KPC-71	1.13	KPC beta-lactamase	1.92	OXA-900	1.44	KPC-71	1.05
OXA-12	0.88	KPC-71	1.13	KPC-71	1.19	OXA-12	0.48
OXA-732	0.66	OXA-780	0.62	OXA-12	0.57	OXA-732	0.41
OXA-46	0.55	OXA-900	0.47	OXA-732	0.20	OXA-780	0.25
OXA-296	0.32	OXA-732	0.30	OXA-780	0.17	OXA-900	0.25
OXA-780	0.32	OXA-46	0.14	OXA-46	0.10	OXA-347	0.20
OXA-347	0.24	OXA-164	0.11	OXA-347	0.09	OXA-46	0.13
OXA-164	0.14	OXA-347	0.11	OXA-164	0.06	OXA-209	0.06
OXA-900	0.08	OXA-372	0.04	OXA-372	0.03	OXA-372	0.03
OXA-209	0.07	OXA-209	0.03	OXA-209	0.02	OXA-164	0.03
OXA-372	0.03	OXA-308	0.01	OXA-308	0.01	OXA-296	0.02
OXA-308	0.01	OXA-464	0.01	OXA-299	0.01	OXA-308	0.01
OXA-664	0.00810	OXA-9	0.01	OXA-664	3.44E-03	OXA-664	0.01
OXA-9	0.00697	OXA-830	0.01	OXA-9	1.98E-03	OXA-9	0.01
OXA-37	0.00389	OXA-664	0.0045	IMP-35	1.94E-03	OXA-60	0.0044
IMP-35	0.00292	OXA-299	0.0043	OXA-830	1.76E-03	OXA-37	0.0034
OXA-18	0.00292	IMP-35	0.0027	IMP beta-lactamase	1.35E-03	OXA-299	0.0029
OXA-830	0.00130	OXA-18	0.0025	OXA-37	8.59E-04	OXA-830	0.0027
OXA-299	0.00097	IMP beta-lactamase	0.0021	OXA-464	7.85E-04	OXA-18	0.0022
OXA-464	0.00097	OXA-37	0.0013	NDM beta-lactamase	5.60E-04	OXA-464	0.0021
IMP beta-lactamase	0.00065	OXA-60	0.0011	OXA-60	5.60E-04	IMP-35	0.0009
OXA-443	0.00065	NDM beta-lactamase	0.0003	OXA-18	5.23E-04	OXA-198	0.0009
OXA-902	0.00065	OXA-198	0.0003	OXA-568	5.23E-04	OXA-568	0.0009
NDM beta-lactamase	0.00049	OXA-902	0.0003	IMP-27	4.86E-04	OXA-85	0.0007
OXA-198	0.00032	OXA-62	0.0003	OXA-198	2.99E-04	OXA-902	0.0005
OXA-154	0.00016	OXA-85	0.0003	OXA-443	2.24E-04	NDM beta-lactamase	0.0003
OXA-243	0.00016	OXA-296	0.0003	OXA-902	1.87E-04	IMP beta-lactamase	0.0003
OXA-60	0.00016	OXA-443	0.0002	OXA-153	1.49E-04	OXA-55	0.0003
		OXA-568	0.0002	OXA-727	1.49E-04	OXA-727	0.0003
		IMP-27	0.0001	OXA-154	1.12E-04		
		OXA-243	0.0001	OXA-243	7.47E-05		
				OXA-45	7.47E-05		
				OXA-258	3.74E-05		
				OXA-85	3.74E-05		

Results: AMR Metagenomic Diversity




WWTP 2



KPC ✓
OXA-48 ✓
IMP ✓
NDM ✗
VIM ✓

WWTP 2							
Carbapenem Resistance Genes of Interest [Abundance (%) Rank]							
Week 1		Week 5		Week 9		Week 13	
OXA beta-lactamase	6.32	OXA beta-lactamase	6.27	OXA beta-lactamase	7.80	OXA beta-lactamase	8.27
KPC beta-lactamase	2.77	KPC beta-lactamase	4.05	KPC beta-lactamase	2.04	KPC beta-lactamase	2.54
KPC-71	1.63	KPC-71	2.44	KPC-71	1.22	KPC-71	1.47
OXA-732	1.10	OXA-12	1.66	OXA-732	0.78	OXA-12	0.72
OXA-12	0.67	OXA-780	0.71	OXA-12	0.66	OXA-732	0.71
OXA-164	0.63	OXA-732	0.50	OXA-347	0.51	OXA-780	0.65
OXA-780	0.52	OXA-347	0.40	OXA-46	0.51	OXA-347	0.54
OXA-347	0.29	OXA-164	0.33	OXA-780	0.45	OXA-46	0.48
OXA-46	0.28	OXA-46	0.31	OXA-164	0.21	OXA-164	0.21
OXA-209	0.09	OXA-209	0.13	OXA-372	0.18	OXA-209	0.17
OXA-372	0.07	OXA-372	0.05	OXA-209	0.15	OXA-372	0.09
OXA-9	0.01	OXA-85	0.03	OXA-9	0.04	OXA-372	0.09
OXA-85	0.01	OXA-9	0.01	OXA-85	0.03	OXA-9	0.02
OXA-664	0.01	OXA-830	0.01	OXA-664	0.01	OXA-85	0.01
OXA-464	0.0043	OXA-902	0.01	OXA-464	0.0034	OXA-308	0.01
OXA-308	0.0037	OXA-464	0.0043	OXA-18	0.0027	OXA-664	0.01
OXA-60	0.0036	OXA-664	0.0024	OXA-902	0.0025	OXA-60	0.01
OXA-443	0.0033	OXA-37	0.0023	OXA-60	0.0022	OXA-443	0.0028
OXA-902	0.0030	OXA-153	0.0011	OXA-308	0.0018	OXA-18	0.0020
OXA-830	0.0027	IMP-35	0.0011	OXA-153	0.0016	OXA-830	0.0019
OXA-37	0.0017	OXA-60	0.0011	OXA-830	0.0015	IMP beta-lactamase	0.0017
IMP-35	0.0009	OXA-443	0.0010	OXA-299	0.0009	OXA-902	0.0017
IMP beta-lactamase	0.0006	IMP beta-lactamase	0.0006	OXA-62	0.0009	OXA-464	0.0016
OXA-153	0.0004	OXA-18	0.0006	IMP-35	0.0007	IMP-35	0.0015
OXA-299	0.0004	OXA-154	0.0005	IMP beta-lactamase	0.0007	OXA-299	0.0015
VIM-7	0.0003	OXA-308	0.0005	OXA-443	0.0007	OXA-62	0.0012
OXA-243	0.0003	OXA-62	0.0005	OXA-45	0.0006	OXA-37	0.0011
OXA-62	0.0003	OXA-299	0.0004	VIM beta-lactamase	0.0003	OXA-568	0.0009
OXA-18	0.0002	OXA-45	0.0003	OXA-900	0.0003	OXA-45	0.0009
VIM beta-lactamase	0.0001	IMP-27	0.0002	OXA-198	0.0003	VIM beta-lactamase	0.0005
OXA-727	0.0001	OXA-215	0.0002	OXA-37	0.0001	OXA-900	0.0003
		OXA-258	0.0002	OXA-243	0.0001	OXA-258	0.0003
		OXA-727	0.0002			OXA-153	0.0001
		VIM-7	0.0001			OXA-243	0.0001
		VIM beta-lactamase	0.0001			OXA-154	0.0001
		OXA-568	0.0001			OXA-198	0.0001
		OXA-243	0.0001				

Results: AMR Metagenomic Diversity




- Both WWTPs showed **very similar AMR genetic diversity** despite being located on **two different islands**.

#IDLabCon



Results: Bacterial Diversity





- Both WWTPs showed that **highly comparable bacterial diversity** despite being located on **two different islands**.
- PCA analysis showed that, **by WWTPs**, bacterial diversity were not distinctly different.

#IDLabCon



Results: Bacterial Diversity



- PCA analysis showed that, **by Weeks**, bacterial diversity were not distinctly different too.
- Especially Week 5

Summary & Conclusion




#IDLabCon

Summary



- In both WWTPs across the State of Hawaii, **KPC** and **OXA-48** were **always detected at a high levels**; whereas **NDM**, **IMP**, and **VIM** were **detected at a low levels or not detected at all**.
- AMR sequencing data were **correlated** with the dPCR data, where the more frequent detected carbapenem resistance genes by dPCR were also showed higher abundances in the sequencing data.
- Bacterial diversity were **not distinctly different** in both WWTPs.

Conclusion

- Both WWTPs showed surprisingly **very similar AMR spatial & temporal distribution, AMR genetic diversity, and bacterial components** despite being located on two different islands.



Takeaways



Human



Environment

Biofilm



Acknowledgement



NATIONAL™
WASTEWATER
SURVEILLANCE
SYSTEM



SPECIAL THANKS

QIAGEN

Bi Linton

Joshua Boston

Shireen Flores, Ph.D.

Leah Riefberg

Sherina Edwards

Tamar Vardi

Matthew Collin, Ph.D.

Shawn Prince

Elizabeth Conzevoy

Ronny Kellner

Dylan Barbera

Ceres Nanosciences

Tara Jones-Roe

Lauren Saunders, Ph.D.

Daniel Goldfarb

HDOH-SLD

Maria Steadmon, Ph.D.

Prakit Saingam, Ph.D.

Carl Min Ki Jeon, Ph.D.

Stella Remigio

Pamela O'Brien

Daniel Strange, Ph.D.


Derek Harauchi

Edward Desmond, Ph.D.

Dayna Ornellas


#IDLabCon





THANK YOU

Questions? Please contact
doris.di.nsw@doh.hawaii.gov



Wastewater Team in HI

